



Tray Layout Guide

Smith+Nephew

MAVERICK 
Unilateral External Fixation




MAVERICK[®] Top Instrument Tray



Reference	Cat No.	Description	Qty
A	71094161	MAVERICK Trocar Short	2
B	71094151	MAVERICK Tissue Protector Short	2
C	71094024	MAVERICK Short Drill 2.4mm for 3mm Pin	2
D	71094315	MAVERICK Half Pin SS 3 x 100 x 15mm	2
E	71094325	MAVERICK Half Pin SS 3 x 100 x 25mm	2
F	71094032	MAVERICK Short Drill 3.2mm for 4mm Pin	2
G	71094420	MAVERICK Half Pin SS 4 x 100 x 20mm	2
H	71094435	MAVERICK Half Pin SS 4 x 100 x 35mm	2
I	71094162	MAVERICK Trocar Medium	2
J	71094152	MAVERICK Tissue Protector Medium	2
K	71094040	MAVERICK Medium Drill 4.0mm for 5mm x 160mm Pin	2
L	71094535	MAVERICK Half Pin SS 5 x 160 x 35mm	4
M	71094555	MAVERICK Half Pin SS 5 x 160 x 55mm	4
N	71094163	MAVERICK Trocar Long	2
O	71094153	MAVERICK Tissue Protector Long	2
P	71094149	MAVERICK Long Drill 4.0mm for 5mm x 200mm & 250mm Pin	2
Q	71094545	MAVERICK Half Pin SS 5 x 200 x 45mm	4
R	71094565	MAVERICK Half Pin SS 5 x 200 x 65mm	4
S	71094585	MAVERICK Half Pin SS 5 x 250 x 85mm	4
T	71094156	MAVERICK Tissue Protector 6mm	2
U	71094166	MAVERICK Trocar 6mm	2
V	71094145	MAVERICK Long Drill 4.5mm for 6mm Pin	2
W	71094655	MAVERICK Half Pin SS 6 x 200 x 55mm	2
X	71094695	MAVERICK Half Pin SS 6 x 250 x 95mm	2

Reference	Cat No.	Description	Qty
Y	71094575	MAVERICK Transfixing Pin SS 5mm x 275mm	2
Z	71094155	MAVERICK Tissue Protector Transfixing Pin	2
ZZ	71094165	MAVERICK Trocar Transfixing Pin	2
AA	71094102	MAVERICK Ratchet Wrench 10mm	1
BB	71094101	MAVERICK AO T-Handle 10mm	1
CC	71094105	MAVERICK T-Handle 10mm Long	1
DD	71094104	MAVERICK Afterburner	1
EE	71094103	MAVERICK Tissue Protector Handle	1

MAVERICK[®] Bottom Instrument Tray



Reference	Cat No.	Description	Qty
A	71094030	MAVERICK Post 30	4
B	71094060	MAVERICK Post Straight	2
C	71094090	MAVERICK Post 90	2
D	71094075	MAVERICK Multi-Pin Clamp 75mm	2
E	71094045	MAVERICK Multi-Pin Clamp 45mm	2
F	71094020	MAVERICK Swivel Clamp 11mm to 5-6mm	12
G	71094015	MAVERICK Combi Clamp 11mm to 5-6mm	4
H	71094111	MAVERICK Swivel Attachment	4
I	71094100	MAVERICK Bar 11mm x 100mm	2
J	71094150	MAVERICK Bar 11mm x 150mm	2
K	71094200	MAVERICK Bar 11mm x 200mm	3
L	71094250	MAVERICK Bar 11mm x 250mm	3
M	71094300	MAVERICK Bar 11mm x 300mm	3
N	71094350	MAVERICK Bar 11mm x 350mm	3
O	71094400	MAVERICK Bar 11mm x 400mm	3
P	71094500	MAVERICK Bar 11mm x 500mm	3
Q	71094210	MAVERICK Curved Bar 11mm x 210mm	1
R	71094380	MAVERICK Curved Bar 11mm x 380mm	1
S	71094190	MAVERICK Bent Bar 11mm x 190mm	1

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, Tennessee 38116
USA

www.smith-nephew.com

°Trademark of Smith+Nephew
All Trademarks acknowledged
©2024 Smith+Nephew
42735 V1 04/24