

Introducing the new **RENASYS[◇] EDGE** Negative Pressure Wound Therapy device

Easy to use:

- Pump offers an intuitive, user-friendly interface for easy to learn operation and troubleshooting¹

Built for efficiency:

- Streamlined cleaning process between patients takes less than 5 minutes²

Patient-friendly:

- Size and weight of device is designed to allow patients to continue with their daily lives and activities³

Smith+Nephew

RENASYS[◇] EDGE

Negative Pressure Wound
Therapy System



The RENASYS[◇] EDGE Device



RENASYS[◇] EDGE
Negative Pressure Wound
Therapy System

Comparing the devices:

The **RENASYS[®] EDGE** pump:

- Smaller footprint*⁴
- More easily portable*⁴
- Has a longer battery life*³

RENASYS TOUCH



RENASYS EDGE



*Compared to RENASYS TOUCH System

	TOUCH	EDGE
Weight	1.22Kg	1.13 Kg
Dimensions	18cmx19cm	15cmx11cm
Footprint area	217 cm ²	165 cm ²
Need annual service and maintenance	Yes	No
Canister with solidifier	300ml or 800ml	300ml hidden or 800ml
Integrated carry handle	No	Yes
Self-standing	Kickstand with 300ml canister	Yes
Display	Full-colour screen	Full-color screen
Operating time	10-16 hours (therapy)	> 24 hours (therapy)
Pressure range	-25 to -200mmHg	-25 to -200mmHg
Therapy settings	Continuous and intermittent	Continuous and Intermittent
Alarm volumes	Selectable High: 74dB Medium: 68dB Low: 60dB	Selectable High: 75dB Medium: 68dB Low: 58dB
Adjustable draw down rates while therapy starts	Yes	Yes
On-board tutorials	Yes	Yes
Pump self-test functionality	No	Yes
O-ring	Serviceable part. Integrated in the pump.	Non-serviceable part. Integrated in disposable canister
Odour filter	Serviceable part. Integrated in the pump.	Non-serviceable part. Integrated in disposable canister
Home healthcare lock	Patient mode restricts access to therapy settings	Pump locks automatically
Autostart after therapy paused for too long	No	Yes – after 1 hour of inactivity
Accessories	Carry bag, carry straps and IV Pole clamp	Carry bag and carry straps
Pump NFC sensor linked to on-line resources	No	Yes

RENASYS[◇] EDGE System ordering codes:

Item description	Code	Carton
RENASYS EDGE device	66803126	1
RENASYS EDGE 300ml canister with solidifier	66803139	5
RENASYS EDGE 300ml canister without solidifier	66803141	5
RENASYS EDGE 800ml canister with solidifier	66803140	5
RENASYS EDGE 800ml canister without solidifier	66803142	5

Item description	Code	Carton
RENASYS EDGE carry bag	66803136	1
RENASYS EDGE carry straps	66803137	1
RENASYS EDGE transit case	66803138	1
RENASYS Y-Connector	66800971	10

Smith & Nephew Pty Ltd
Australia
T +61 2 9857 3999
F +61 2 9857 3900
smith-nephew.com/en-au

Smith & Nephew Ltd
New Zealand
T +64 9 820 2840
F +64 9 820 2841
smith-nephew.com/en-nz

◇Trademark of Smith+Nephew
All trademarks acknowledged
© June 2024 Smith+Nephew
41043-anz V1 06/24

This material is intended for healthcare professionals. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Shaping 
**what's possible
in wound care**

References: **1.** Smith+Nephew 2022. RENASYS EDGE System Human Factors Summative Report Summary. Internal Report. CSD.AWM.22.071. **2.** Smith+Nephew 2022. Summary of RENASYS EDGE pump cleaning, self-test and maintenance. Internal Report. CSD.AWM.22.068. **3.** Smith+Nephew 2022. How the RENASYS EDGE Negative Pressure Wound Therapy System provides continuity of care to the patient Internal Report. EO.AWM.PCS270.003.V1. **4.** Smith+Nephew 2022. Summary of footprint, portability, wearability, weight and audible noise for the RENASYS[◇] EDGE system. Internal Report. CSD.AWM.22.067.