Simply advanced

The EVOS WRIST System provides a comprehensive solution for treating simple to complex wrist fractures.


Smith&nephew EVOS* WRIST Plating System

Supporting healthcare professionals

Wrist fractures are complicated enough. Why be limited by your system?

EVOS* WRIST Plating System offers:


- Flexibility to choose between stainless and titanium volar options
- Plates and screws designed to ensure low-profile constructs whether screws are on or up to 15° off-axis
- Plating options that considers all distal radius and distal ulna fracture needs
- Titanium dorsal and Stainless steel forearm options



Expansive plating system







Seamless compatibility

EVOS° WRIST Titanium and VLP Mini Mod have you covered from forearm to finger tip. – Compatible screws, drills, and drivers between systems.



EVOS° WRIST Stainless and EVOS° MINI offer stainless options for forearm, wrist, and hand fixation. – Compatible screws, drills, and drivers between systems.



Advanced technology to handle simple to complex wrist fractures

Titanium Screws



Stainless Steel Screws



Variable-angle locking technology





One drill bit One driver

Stainless

1.8mm Drill T7 Driver





1.8mm Pegs

2.5mm Drill2.5mm Hex Driver



3.5mm Screws



1.8mm Drill T7 Driver

Titanium



2.4mm Screws


1.8mm Pegs



3.5mm Screws


Low-profile constructs

Plates and screws designed to ensure low profile constructs whether screws are on or up to 15° off-axis.




Evolutionary approach to plate and screw designs

Anatomic Fracture Specific Dorsal Options




Dorsal Buttress Plate Standard

Dorsal Buttress Plate Wide



Volar options for intra-articular and extra-articular fractures






Anatomical Forearm Plates


Proximal Radial Shaft Plate

Radial Shaft Plate




EVOS^{*} WRIST Titanium System

All inclusive designs that include plates, screws & instrumentation needed for volar plating cases.







Note: Works with EVOS WRIST Titanium system.



EVOS WRIST Stainless System







Note: Works with EVOS SMALL System



EVOS SMALL Instrument Tray



EVOS SMALL Screw Tray

Smith & Nephew Pty Ltd Australia T +61 2 9857 3999 F +61 2 9857 3900 Smith & Nephew Ltd New Zealand T +64 9 820 2840 F +64 9 820 2841 www.smith-nephew.com/australia www.smith-nephew.com/new-zealand °Trademark of Smith+Nephew. All trademarks acknowledged. 19085-anz V1 02/21 For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.