



Versatility for tonsillectomy and adenoidectomy

 **smith&nephew**
COBLATION[◇]
PROCISE[◇]
Adenotonsillectomy Wands

PROCISE[◇] XP and PROCISE EZ adenotonsillectomy Wands

Key features

Robust active electrode configuration

- Removes tissue for both tonsillectomy and adenoidectomy
- The PROCISE XP Wand tip features triple-wire active electrodes
- The PROCISE EZ Wand tip features a flat screen active electrode

Decreased diameter on working shaft of Wand

- Increased surgical field visualization – the PROCISE EZ Wand tip features the smallest shaft diameter among all ENT COBLATION[®] Wands
- Same working length as EVAC[°] Wands
- Bendable shafts for easy access to choanae during adenoidectomy

Unique saline delivery system

- Provides tip with saline regardless of orientation
- Unobstructed port clears surgical field to reduce clogging



Ordering information

PROCISE Wands

Reference #	Description
EICA8872-01	PROCISE XP WAND
EICA8870-01	PROCISE EZ WAND
EC8000-01	COBLATOR [®] II CONTROLLER

ArthroCare Corporation

7000 West William Cannon Drive
Austin, TX 78735
USA

www.smith-nephew.com

Order Entry: 1-800-343-5717
Order Entry Fax: 1-888-994-2782

© 2015 Smith & Nephew, Inc.
*Trademark of Smith & Nephew. Reg. US Pat. & TM Office.

