



Control the **flow** for optimal performance
on multiple tissue types

 **smith&nephew**

COBLATION[◇]


FLOW 50[◇]

Wand with FLOW~IQ[◇]
Technology



Control the **flow** for optimal performance on multiple tissue types

- FLOW~IQ[®] Technology provides automated outflow control that delivers a new standard of ablation performance.
- The FLOW 50[°] Wand adapts to multiple tissue types such as articular cartilage, meniscus, and ligamentous tissue, efficiently resecting targeted tissue while preserving healthy surrounding tissue.



Cartilage treatment



Meniscal sculpting



ACL debridement

Technical specifications

Catalog number	72290037
Controller	RF20000a or RF20000 Systems only
Shaft Diameter	3.76 mm (0.148")
Working length	135 mm (5.3")
Tip Width	2.95 mm (0.116")
Default set point	Med



Improved Integrated Finger Switch (IFS) Technology

Ordering information

Reference #	Description
72290037	WEREWOLF [®] FLOW 50 COBLATION [®] Wand



Slim electrode design facilitates access to targeted tissue

Supporting healthcare professionals for over 150 years

ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735
USA

www.smith-nephew.com
Order Entry: 1-800-343-5717
Order Entry Fax: 1-888-994-2782