# An absorbent<sup>1</sup> dressing for the management of light-to-moderately exuding wounds

MELOLIN<sup>•</sup> Dressing helps protect the wound from physical contaminants.<sup>2</sup>



MELOLIN<sup>6</sup> Low Adherent Absorbent Wound Dressing consists of an absorbent<sup>1</sup> cotton and polyester fibre pad which is heat bonded to a thin low adherent perforated polyester film.

# Features and benefits

# Designed to be low adherent and absorbent MELOLIN Dressing is designed to be low

adherent and absorbent.<sup>1</sup>

# Protects the wound<sup>2</sup>

MELOLIN Dressing has a long history of medical use<sup>3,4</sup> and protects the wound from physical contaminants.<sup>2</sup>

### Water repellent backing<sup>5</sup>

MELOLIN Dressing has a water repellent backing, which minimises fluid strikethrough whilst dressing is in place.<sup>5</sup>

# Non-sensitising<sup>6</sup>\* and breathable<sup>7</sup>

MELOLIN Dressing is non-sensitising<sup>6\*</sup> and allows the skin to breathe.<sup>7</sup>

# Remains integrity when cut<sup>8,9</sup>

MELOLIN Dressing is easy to cut to shape and retains its integrity when cut.8,9

# Compatibility

MELOLIN Dressing can be used in conjunction with OPSITE<sup>4</sup> FLEXIFIX, OPSITE FLEXIFIX Gentle and OPSITE FLEXIGRID.9,10

## Indications

For the management of a wide variety of light-to-moderately exuding wounds including clean sutured wounds, abrasions, lacerations and minor burns.



#### **MELOLIN** Dressings

Non sterile				Sterile		
S+N Code	Size	Carton		S+N Code	Size	Carton
66974942	10cm x 10cm	150		66974932	5cm x 5cm	25
66974943	10cm x 20cm	75		66974940	5cm x 5cm	100
66974944	20cm x 30cm	25		66974933	10cm x 10cm	10
66974930	50cm x 7m	1 roll		66974941	10cm x 10cm	100
				66974939	10cm x 20cm	100

#### Advanced Wound Management www.smith-nephew.com

♦Trademark of Smith+Nephew All Trademarks acknowledged ©May 2020 Smith+Nephew AWM-AWD-17616. GMC1063

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

#### \*As demonstrated in vivo

T +44 (0) 1482 225181

F+44 (0) 1482 328326

Smith & Nephew Medical Ltd

101 Hessle Road Hull HU3 2BN

#### References

UK

 Smith & Nephew 2008. Generation of claims support for MELOLIN. Internal report. CCA 1055/R1.
Smith & Nephew 2019. Protects the wound. Internal report. EO.AWM.PCS0101.001.v1.
Robb W. Clinical trial of MELOLIN: A new non-adherent dressing. British Journal of Plastic Surgery. 1961;14(1):47-49.
Gentle M. MELOLIN trial. Nursing Mirror. 1970:1-4.
Smith & Nephew 2013. Safety summary. Internal report. 7. Smith & Nephew 2019. Moisture vapour permeability. testing for MELOLIN / CUTILIN<sup>®</sup>. Internal report. U/071/R1 V1. **8**. Smith & Nephew 2019. Fibre shedding testing for MELOLIN - CUTILIN. Internal report. U/071/R3 V2. **9**. Smith & Nephew 2019. Verification testing of usability related design input for MELOLIN / CUTILIN. Internal report. U/071/R4 V2. **10**. Smith & Nephew 2019. Physical compatibility testing for MELOLIN. Internal report. U/071/R6.