

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 1. Identification

#### 1.1. Product identifier

**Product Identity** No Sting Skin Prep: Wipes and Swabs

**Alternate Names** No Sting Skin Prep: Wipes and Swabs

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

**Intended use** Skin Preparation

**Application Method** See Technical Data Sheet.

#### 1.3. Details of the supplier of the safety data sheet

**Company Name** Smith & Nephew  
970 Lake Carillon Drive, Suite 110  
St. Petersburg, FL 33716

#### Emergency

**Customer Service: Smith & Nephew** 1-800-876-1261

### 2. Hazard(s) identification

#### 2.1. Classification of the substance or mixture

Flam. Liq. 2;H225 Highly Flammable liquid and vapor.

Aquatic Acute 1;H400 Very toxic to aquatic life.

Aquatic Chronic 2;H411 Toxic to aquatic life with long lasting effects.

#### 2.2. Label elements

Using the Toxicity Data listed in section 11 and 12 the product is labeled as follows.



**Danger**

H225 Highly flammable liquid and vapor.

H400 Very toxic to aquatic life.

H411 Toxic to aquatic life with long lasting effects.

#### [Prevention]:

P210 Keep away from heat / sparks / open flames / hot surfaces - No smoking.

P235 Keep cool.

P240 Ground / bond container and receiving equipment.

P241 Use explosion-proof electrical / ventilating / light / equipment.

P242 Use only non-sparking tools.

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

P243 Take precautionary measures against static discharge.

P273 Avoid release to the environment.

P280 Wear protective gloves / eye protection / face protection.

**[Response]:**

P303+361+353 IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower.

P370+378 In case of fire: Use extinguishing media listed in section 5 of SDS for extinction.

P391 Collect spillage.

**[Storage]:**

P403+233 Store in a well ventilated place. Keep container tightly closed.

**[Disposal]:**

P501 Dispose of contents / container in accordance with local / national regulations.

### 3. Composition/information on ingredients

This product contains the following substances that present a hazard within the meaning of the relevant State and Federal Hazardous Substances regulations.

Ingredient/Chemical Designations	Weight %	GHS Classification	Notes
Disiloxane, hexamethyl- CAS Number: 0000107-46-0	75 - 100	Flam. Liq. 2;H225 Aquatic Acute 1;H400 Aquatic Chronic 2;H411	[1]
Acrylate Copolymer CAS Number: Proprietary	10 - 25	Not Classified	[1]

In accordance with paragraph (i) of §1910.1200, the specific chemical identity and/or exact percentage (concentration) of composition has been withheld as a trade secret.

[1] Substance classified with a health or environmental hazard.

[2] Substance with a workplace exposure limit.

[3] PBT-substance or vPvB-substance.

\*The full texts of the phrases are shown in Section 16.

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 4. First aid measures

#### 4.1. Description of first aid measures

<b>General</b>	In all cases of doubt, or when symptoms persist, seek medical attention. Never give anything by mouth to an unconscious person.
<b>Inhalation</b>	Remove to fresh air, keep patient warm and at rest. If breathing is irregular or stopped, give artificial respiration. If unconscious place in the recovery position and obtain immediate medical attention. Give nothing by mouth.
<b>Eyes</b>	Irrigate copiously with clean water for at least 15 minutes, holding the eyelids apart and seek medical attention.
<b>Skin</b>	Remove contaminated clothing. Wash skin thoroughly with soap and water or use a recognized skin cleanser.
<b>Ingestion</b>	If swallowed obtain immediate medical attention. Keep at rest. Do NOT induce vomiting.

#### 4.2. Most important symptoms and effects, both acute and delayed

<b>Overview</b>	Flammable. May cause eye, skin and respiratory tract irritation. See section 2 for further details.
-----------------	---

### 5. Fire-fighting measures

#### 5.1. Extinguishing media

Extinguishing powder, sand, alcohol-resistant foam, carbon dioxide

#### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition: Oxides of Carbon

Keep away from heat / sparks / open flames / hot surfaces - No smoking.

Keep cool.

Ground / bond container and receiving equipment.

Use explosion-proof electrical / ventilating / light / equipment.

Use only non-sparking tools.

Take precautionary measures against static discharge.

#### 5.3. Advice for fire-fighters

Use self-contained breathing apparatus and wear protective clothing. Cool endangered containers with water spray jet. Suppress gases/vapors/mists with water spray jet.

**ERG Guide No.** 133

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 6. Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

Put on appropriate personal protective equipment (see section 8).

#### 6.2. Environmental precautions

Do not allow spills to enter drains or waterways.

Use good personal hygiene practices. Wash hands before eating, drinking, smoking or using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

#### 6.3. Methods and material for containment and cleaning up

Isolate area and keep unnecessary personnel away.

Pick up with absorbent material (e.g. sand, acid binder, universal binder, sawdust). Send in suitable containers for recover or disposal.

### 7. Handling and storage

#### 7.1. Precautions for safe handling

Vapors can form an explosive mixture with air. Take precautionary measures against static charges. Keep away from sources of heat and ignition. Use explosion-proof equipment/fittings and non-sparking tools.

See section 2 for further details. - [Prevention]:

#### 7.2. Conditions for safe storage, including any incompatibilities

Handle containers carefully to prevent damage and spillage.

Incompatible materials: No data available.

Containers which are opened must be carefully resealed and kept upright to prevent leakage. Always keep in containers of same material as the original one. Keep container tightly closed and dry in a cool, well-ventilated place. Protect from heat and direct sunlight.

See section 2 for further details. - [Storage]:

#### 7.3. Specific end use(s)

No data available.

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 8. Exposure controls and personal protection

#### 8.1. Control parameters

##### Exposure

CAS No.	Ingredient	Source	Value
0000107-46-0	Disiloxane, hexamethyl-	OSHA	No Established Limit
		ACGIH	No Established Limit
		NIOSH	No Established Limit
		Supplier	No Established Limit
Proprietary	Acrylate Copolymer	OSHA	No Established Limit
		ACGIH	No Established Limit
		NIOSH	No Established Limit
		Supplier	No Established Limit

#### Carcinogen Data

CAS No.	Ingredient	Source	Value
0000107-46-0	Disiloxane, hexamethyl-	OSHA	Select Carcinogen: No
		NTP	Known: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;
Proprietary	Acrylate Copolymer	OSHA	Select Carcinogen: No
		NTP	Known: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;

#### 8.2. Exposure controls

##### Respiratory

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of aerosol and mist formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

##### Eyes

Wear safety glasses with side shields.

##### Skin

Wear appropriate clothing for the industrial setting.

##### Engineering Controls

Provide adequate ventilation. Where reasonably practicable this should be achieved by the use of local exhaust ventilation and good general extraction. If these are not sufficient to maintain concentrations of particulates and any vapor below occupational exposure limits suitable respiratory protection must be worn.

##### Other Work Practices

Use good personal hygiene practices. Wash hands before eating, drinking, smoking or using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

See section 2 for further details. - [Prevention]:

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 9. Physical and chemical properties

<b>Appearance</b>	Colorless Liquid
<b>Odor</b>	Odorless
<b>Odor threshold</b>	Not determined
<b>pH</b>	Not determined
<b>Melting point / freezing point</b>	Not applicable
<b>Initial boiling point and boiling range</b>	37°C
<b>Flash Point</b>	-1°C
<b>Evaporation rate (Ether = 1)</b>	Not determined
<b>Flammability (solid, gas)</b>	Not Applicable
<b>Upper/lower flammability or explosive limits</b>	<b>Lower Explosive Limit:</b> Not determined <b>Upper Explosive Limit:</b> Not determined
<b>Vapor pressure (Pa)</b>	Not determined
<b>Vapor Density</b>	0.76 g/cm <sup>3</sup>
<b>Specific Gravity</b>	Not determined
<b>Solubility in Water</b>	Insoluble
<b>Partition coefficient n-octanol/water (Log Kow)</b>	Not Measured
<b>Auto-ignition temperature</b>	Not determined
<b>Decomposition temperature</b>	Not available
<b>Viscosity (cSt)</b>	Not available
<b>VOC Content</b>	Not determined

#### 9.2. Other information

No other relevant information.

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 10. Stability and reactivity

#### 10.1. Reactivity

Hazardous Polymerization will not occur.

#### 10.2. Chemical stability

Stable under normal circumstances.

#### 10.3. Possibility of hazardous reactions

No data available.

#### 10.4. Conditions to avoid

No data available.

#### 10.5. Incompatible materials

No data available.

#### 10.6. Hazardous decomposition products

Oxides of Carbon

### 11. Toxicological information

#### Acute toxicity

Ingredient	Oral LD50, mg/kg	Skin LD50, mg/kg	Inhalation Vapor LC50, mg/L/4hr	Inhalation Dust/Mist LC50, mg/L/4hr	Inhalation Gas LC50, ppm
Disiloxane, hexamethyl- - (107-46-0)	5,000.00, Rat - Category: 5	5,000.00, Rabbit - Category: 5	106.00, Rat - Category: NA	No data available	No data available
Acrylate Copolymer - (Proprietary)	No data available	No data available	No data available	No data available	No data available

Note: When no route specific LD50 data is available for an acute toxin, the converted acute toxicity point estimate was used in the calculation of the product's ATE (Acute Toxicity Estimate).

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

Classification	Category	Hazard Description
Acute toxicity (oral)	---	Not Applicable
Acute toxicity (dermal)	---	Not Applicable
Acute toxicity (inhalation)	---	Not Applicable
Skin corrosion/irritation	---	Not Applicable
Serious eye damage/irritation	---	Not Applicable
Respiratory sensitization	---	Not Applicable
Skin sensitization	---	Not Applicable
Germ cell mutagenicity	---	Not Applicable
Carcinogenicity	---	Not Applicable
Reproductive toxicity	---	Not Applicable
STOT-single exposure	---	Not Applicable
STOT-repeated exposure	---	Not Applicable
Aspiration hazard	---	Not Applicable

## 12. Ecological information

### 12.1. Toxicity

Very toxic to aquatic life.

Toxic to aquatic life with long lasting effects.

No additional information provided for this product. See Section 3 for chemical specific data.

### Aquatic Ecotoxicity

Ingredient	96 hr LC50 fish, mg/l	48 hr EC50 crustacea, mg/l	ErC50 algae, mg/l
Disiloxane, hexamethyl- - (107-46-0)	3.02, Oncorhynchus mykiss	314.00, Daphnia magna	Not Available
Acrylate Copolymer - (Proprietary)	Not Available	Not Available	Not Available

### 12.2. Persistence and degradability

There is no data available on the preparation itself.

### 12.3. Bioaccumulative potential

Not Measured

### 12.4. Mobility in soil

No data available.

### 12.5. Results of PBT and vPvB assessment

This product contains no PBT/vPvB chemicals.

### 12.6. Other adverse effects

No data available.

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 13. Disposal considerations

**13.1. Waste treatment methods**

Observe all federal, state and local regulations when disposing of this substance.

### 14. Transport information

**Notes:**

1. May qualify for transportation under 'Limited Quantities' and/or 'Excepted Quantities' provisions. Consult applicable DOT, IMDG, IATA/ICAO regulations depending on transportation mode.
2. May qualify for transportation under 'Consumer Commodity ORM-D' provisions. Consult applicable DOT, IMDG, IATA/ICAO regulations depending on transportation mode.

	<b>DOT (Domestic Surface Transportation)</b>	<b>IMO / IMDG (Ocean Transportation)</b>	<b>ICAO/IATA</b>
<b>14.1. UN number</b>	UN3175	UN3175	UN3175
<b>14.2. UN proper shipping name</b>	UN3175, Solids containing flammable liquid, n.o.s., (Hexamethyldisiloxane), 4.1, II	Solids containing flammable liquid, n.o.s., (Hexamethyldisiloxane)	Solids containing flammable liquid, n.o.s., (Hexamethyldisiloxane)
<b>14.3. Transport hazard class(es)</b>	<b>DOT Hazard Class:</b> 4.1	<b>IMDG:</b> 4.1 <b>Sub Class:</b> Not Applicable	<b>Air Class:</b> 4.1
<b>14.4. Packing group</b>	II	II	II
<b>14.5. Environmental hazards</b>	IMDG Marine Pollutant: Yes ( Disiloxane, hexamethyl- )		
<b>14.6. Special precautions for user</b>	No further information		

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 15. Regulatory information

**Regulatory Overview** The regulatory data in Section 15 is not intended to be all-inclusive, only selected regulations are represented.

**Toxic Substance Control Act ( TSCA)** All components of this material are either listed or exempt from listing on the TSCA Inventory.

**WHMIS Classification** B2

**US EPA Tier II Hazards**

**Fire:** Yes

**Sudden Release of Pressure:** No

**Reactive:** No

**Immediate (Acute):** No

**Delayed (Chronic):** No

**EPCRA 311/312 Chemicals and RQs:**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**EPCRA 302 Extremely Hazardous:**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**EPCRA 313 Toxic Chemicals:**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**Proposition 65 - Carcinogens (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**Proposition 65 - Developmental Toxins (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**Proposition 65 - Female Repro Toxins (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**Proposition 65 - Male Repro Toxins (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**New Jersey RTK Substances (>1%) :**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

**Pennsylvania RTK Substances (>1%) :**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

# Safety Data Sheet

## No Sting Skin Prep: Wipes and Swabs

SDS Revision Date:

06/03/2015

### 16. Other information

The information and recommendations contained herein are based upon data believed to be correct. However, no guarantee or warranty of any kind, expressed or implied, is made with respect to the information contained herein. We accept no responsibility and disclaim all liability for any harmful effects which may be caused by exposure to our products. Customers/users of this product must comply with all applicable health and safety laws, regulations, and orders.

The full text of the phrases appearing in section 3 is:

H225 Highly flammable liquid and vapor.

H400 Very toxic to aquatic life.

H411 Toxic to aquatic life with long lasting effects.

Not Classified

**This is the first version in the GHS SDS format. Listings of changes from previous versions in other formats are not applicable.**

End of Document