

+ ALLEVYN[◇] TRACHEOSTOMY Foam Dressing

For moderate to highly exuding wounds



Smith+Nephew

ALLEVYN[◇]
TRACHEOSTOMY
Foam Dressing

Helping you get **CLOSER TO ZERO[◇]**
human and economic consequence of wounds
smith-nephew.com/allevyn

+ ALLEVYN[◇] TRACHEOSTOMY Foam Dressing

is a soft non-adherent dressing with a 'key hole' aperture to allow use with a tracheostomy tube or other drain or stoma.^{1,2}

Features and benefits

Non-adherent

ALLEVYN TRACHEOSTOMY Dressing is non-adherent to help minimise trauma to the wound and surrounding skin during dressing changes.²

Absorbent

ALLEVYN TRACHEOSTOMY Dressing absorbs exudate to help minimise the risk of maceration to healthy skin.^{3,4}

Prevents leakage

ALLEVYN TRACHEOSTOMY Dressing has a waterproof outer film which assists in preventing leakage of secretions through the front of the dressing.^{4,5}

Easy to use

ALLEVYN TRACHEOSTOMY Dressing is easy to apply, remove and can be cut to size.^{2,3,6}

Bacterial barrier

ALLEVYN TRACHEOSTOMY Dressing has a bacterial barrier top film that aids in the prevention of bacterial contamination.⁷

*As demonstrated *in vitro*.

Indications

Absorbent apertured dressing for the stoma created by tracheostomy or for wound drainage sites.

Do not use ALLEVYN TRACHEOSTOMY Dressings with oxidising agents such as hypochlorite solutions (eg. EUSOL) or hydrogen peroxide, as these can break down the absorbent polyurethane component of the dressing. The scale on the pouch is for guidance to assist dressing selection only, and care should be taken to avoid cross contamination.



ALLEVYN TRACHEOSTOMY Dressings order references

Sizes	Code	Carton
9cm x 9cm	66800740	10

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Advanced Wound Management
Smith & Nephew Medical Ltd
Building 5, Hatters Lane - Watford
Hertfordshire WD18 8YE - UK

www.smith-nephew.com
T +44 (0) 1923 477100
F +44 (0) 1923 477101

◇ Trademark of Smith+Nephew
All Trademarks acknowledged
© July 2020 Smith+Nephew
AWM-AWC-25857 | GMC1128

References

1. Smith & Nephew 2019. Volunteer evaluation of AWM product features. Internal report EO.AWM.PCSgen.007.v1. 2. Smith & Nephew 2006. A prospective user evaluation on the effect of process changes to the manufacture of ALLEVYN NON-ADHESIVE wound dressings. Internal report CER/CE/022/ANA. 3. Kurring PA, Roberts CD, Quinlan D. Evaluation of a hydrocellular dressing in the management of exuding wounds in the community. *British Journal of Nursing* 1994;3(20):1049-1050, 1052-1043. 4. Smith & Nephew 2009. ALLEVYN NON-ADHESIVE Dressing physical properties. Internal report DS/09/013/R3. 5. Smith & Nephew July 2006. Evaluation of leakage at the edges of sealed and shaped ALLEVYN NON-ADHESIVE Dressings. Internal report. DOF 0603011. 6. Smith & Nephew 2015. Cutting of ALLEVYN variants. Internal report DS/14/318/R. 7. Smith & Nephew 2014. Bacterial barrier testing (Wet-Wet) of ALLEVYN NON-ADHESIVE Dressings against *S. marcescens* over 7 days. Internal report DOF 1401002.