

For precise, reliable control
of your instability repairs



Providing precise anatomic fixation of soft tissues

Precise, reliable tissue shift and anatomic repair is now achievable with the SPEEDLOCK® Knotless Suture Anchor. This knotless anchor provides independent suture-lock and accurate suture tensioning after anchor insertion.

The SPEEDLOCK Anchor features independent Suture Lock, which allows you to make incremental micro-tensioning suture adjustments for a precise repair. The internal suture locking mechanism eliminates the bone and suture lock variability that may arise from differences in bone quality and arthroscopic knot tying.

- **Inline drill guide** system provides precise, ergonomic, reliable drill placement.
- **Strong fixation:** The SPEEDLOCK Knotless Suture Anchor demonstrates pullout strength of 250N¹.
- **Independent suture tensioning** allows the surgeon to precisely and reliably reduce the tissue to the underlying bone after anchor insertion.
- **Internal suture lock** allows for precise, reproducible suture tensioning, independent of anchor and bone interference.
- **PEEK implant** is radiolucent, MRI-safe, and generally revisable.



1. In 20 lb. bone foam as documented in report 25374-01 Rev. G for the SPEEDLOCK Knotless Suture Anchor.

Ordering Information

SPEEDLOCK Knotless Suture Anchor

Reference #	Description
OM-7500	SPEEDLOCK Knotless Suture Anchor

SPEEDLOCK Inline Drill System

Reference #	Description
OM-9816	3.5mm Inline Drill Guide
OM-9817	SPEEDLOCK Obturator for Inline Drill Guide
OM-9850	SPEEDLOCK 3.0mm Tapered Twist Drill
OM-9860	SPEEDLOCK 3.0mm Tapered Twist Drill (hard bone)

SPEEDLOCK Drill Guide

Reference #	Description
OM-9715	SPEEDLOCK Drill Guide, 3.5mm, high visibility
OM-9716	SPEEDLOCK Drill Guide, 3.5mm, low profile
OM-9717	SPEEDLOCK Sharp-tipped Obturator
OM-9750	SPEEDLOCK 3.0mm Tapered Twist Drill



Supporting healthcare professionals for over 150 years

Manufactured by:
ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735 USA

Distributed by:
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810 USA

www.smith-nephew.com
T +1 978 749 1000
US Customer Service:
+1 800 343 5717

*Trademark of Smith & Nephew.
©2017 Smith & Nephew. All rights reserved.
Printed in USA. 05616 V2 PN 74027 Rev. B 05/17