

Smith+Nephew REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)
Issue date: 8/7/2026

SECTION 1 Identification

1.1. Product identifier

Product form : Mixture
Product name : REMOVE Adhesive Remover (bottled)
Product code : 403300, 403379

1.2. Other means of identification

No additional information available

1.3. Recommended use of the chemical and restrictions on use

Intended for general public
Use of the substance/mixture : Medical devices
Recommended use : Intended for using on skin, for universal adhesive removal.

1.4. Supplier's details

Smith & Nephew Medical Ltd
101 Hessle Road
Hull, East Ridings of Yorkshire, HU3 2BN
T +44 (0) 1482 225181
<https://www.smith-nephew.com/en/contact-us#general>

1.5. Emergency phone number

Emergency number : 911 (24 Hours)1-800-424-9300 (Chemtrec 24 Hours)

SECTION 2 Hazard Identification

2.1. Classification of the substance or mixture

GHS US classification

Flammable liquid, Category 3
Aspiration hazard, Category 1
Full text of H statements : see section 16

H304 H226 Flammable liquid and vapor.
May be fatal if swallowed and enters airways.

2.2. Label elements

GHS US labeling

Hazard pictograms (GHS US)



Signal word (GHS US) : Danger
Hazard statements (GHS US) : H226 - Flammable liquid and vapor
H304 - May be fatal if swallowed and enters airways
Precautionary statements (GHS US) : P101 - If medical advice is needed, have product container or label at hand.
P102 - Keep out of reach of children.
P103 - Read label before use.

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233 - Keep container tightly closed.
P240 - Ground/Bond container and receiving equipment.
P241 - Use explosion-proof equipment.
P242 - Use non-sparking tools.
P243 - Take action to prevent static discharges.
P301+P310 - If swallowed: Immediately call a poison center or doctor.
P331 - Do NOT induce vomiting.
P370+P378 - In case of fire: Use appropriate media to extinguish.
P403+P235 - Store in a well-ventilated place. Keep cool.
P405 - Store locked up.
P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

2.3. Hazards associated with known or reasonably anticipated uses

No additional information available

2.4. Hazards not otherwise classified

No additional information available

2.5. Unknown acute toxicity

No additional information available

SECTION 3 Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	GHS US classification
C9-12-ISOALKANES	CAS-No.: 90622-57-4	≥ 20 – < 50	Flam. Liq. 3, H226 Asp. Tox. 1, H304
EUCALYPTUS CITRODORA OIL	CAS-No.: 8000-48-4	< 0.1	Eye Irrit. 2, H319 Aquatic Chronic 2, H411 Skin Irrit. 2, H315 Skin Sens. 1, H317

Full text of hazard classes and H-statements : see section 16

SECTION 4 First aid measures

4.1. Description of necessary first-aid measures

First-aid measures general : If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation : May cause drowsiness or dizziness. If experiencing respiratory symptoms: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell.
First-aid measures after skin contact : If skin irritation or rash occurs: Gently wash with plenty of soap and water.

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

First-aid measures after eye contact	: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
First-aid measures after ingestion	: Rinse mouth, do not induce vomiting. (If vomiting occurs, keep head low so that stomach content doesn't get into lungs.).

4.2. Most important symptoms/effects, acute and delayed

Potential Adverse human health effects and symptoms	: Based on available data, the classification criteria are not met.
Symptoms/effects	: Not expected to present a significant hazard under anticipated conditions of normal use.
Symptoms/effects after inhalation	: May cause respiratory irritation.
Symptoms/effects after skin contact	: May cause slight irritation.
Symptoms/effects after eye contact	: May cause slight irritation. stinging. redness, itching, tears.
Symptoms/effects after ingestion	: Ingestion may cause nausea, vomiting and diarrhea. Risk of lung edema.

4.3. Indication of immediate medical attention and special treatment needed, if necessary

Other medical advice or treatment	: Treat symptomatically.
-----------------------------------	--------------------------

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Alcohol-resistant foam. Dry powder. Foam. Water spray or fog. Carbon dioxide.
Unsuitable extinguishing media	: None.

5.2. Specific hazards arising from the chemical

Fire hazard	: Highly flammable liquid and vapor.
Explosion hazard	: May form flammable/explosive vapor-air mixture.
Hazardous decomposition products in case of fire	: Toxic fumes may be released.

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions	: Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.
Protection during firefighting	: Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6 Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures	: None under normal use.
------------------	--------------------------

For non-emergency personnel

Emergency procedures	: Evacuate unnecessary personnel.
----------------------	-----------------------------------

For emergency responders

Protective equipment	: Equip cleanup crew with proper protection.
Emergency procedures	: Ventilate area.

Environmental precautions	: Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.
---------------------------	---

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

6.2. Methods and materials for containment and cleaning up

- | | |
|-------------------------|---|
| Methods for cleaning up | : Take up liquid spill into absorbent material. |
| Other information | : If spilled, may cause the floor to be slippery. |

See Heading 8. Exposure controls and personal protection. For further information refer to section 13.

SECTION 7 Handling and storage

7.1. Precautions for safe handling

- | | |
|-----------------------------------|--|
| Precautions for safe handling | : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor. No open flames. No smoking. Avoid breathing mist, spray. |
| Hygiene measures | : Wash hands, forearms and face thoroughly after handling. |
| Additional hazards when processed | : Handle empty containers with care because residual vapors are flammable. |

7.2. Conditions for safe storage, including incompatibilities

- | | |
|------------------------|---|
| Technical measures | : Proper grounding procedures to avoid static electricity should be followed. Ground/bond container and receiving equipment. Use explosion-proof electrical/ventilating/lighting equipment. |
| Storage conditions | : Keep only in the original container in a cool, well-ventilated place away from: Keep in fireproof place. Keep container tightly closed. |
| Incompatible products | : Strong bases. Strong acids. |
| Incompatible materials | : Direct sunlight. |
| Storage temperature | : < 25 °C |

SECTION 8 Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Appropriate engineering controls

No additional information available

8.3. Individual protection measures, such as personal protective equipment

Personal protective equipment:

Avoid all unnecessary exposure.

Hand protection:
Not required for normal conditions of use
Eye protection:
Not required for normal conditions of use
Respiratory protection:
Not required for normal conditions of use

Other information:

Do not eat, drink or smoke during use.

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

SECTION 9 Physical and chemical properties

9.1. Basic physical and chemical properties

Physical state	: Liquid
Appearance	: Clear, colorless liquid.
Color	: Colourless liquid.
Odor	: pleasant
Odor threshold	: No data available
pH	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Flammability (solid, gas)	: No data available
Vapor pressure	: No data available
Relative vapor density at 20°C	: No data available
Relative density	: No data available
Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Explosion limits	: No data available
Particle characteristics	: No data available

9.2. Data relevant with regard to physical hazard classes (supplemental)

No additional information available

SECTION 10 Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Not established.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures. Open flame. Overheating. Heat. Sparks.

10.5. Incompatible materials

Strong acids. Strong bases.

10.6. Hazardous decomposition products

Fume. Carbon monoxide. Carbon dioxide. May release flammable gases.

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

SECTION 11 Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral) : No data available
Acute toxicity (dermal) : No data available
Acute toxicity (inhalation) : No data available

EUCALYPTUS CITRODORA OIL (8000-48-4)

LD50 oral	2726 mg/kg body weight
LD50 dermal	2480 mg/kg body weight
ATE US (oral)	2726 mg/kg body weight
ATE US (dermal)	2480 mg/kg body weight

Skin corrosion/irritation : No data available

Serious eye damage/irritation : No data available

Respiratory or skin sensitization : No data available

Germ cell mutagenicity : No data available

Carcinogenicity : No data available

Reproductive toxicity : No data available

STOT-single exposure : No data available

STOT-repeated exposure : No data available

Aspiration hazard : May be fatal if swallowed and enters airways.

Potential Adverse human health effects and symptoms : Based on available data, the classification criteria are not met.

Symptoms/effects : Not expected to present a significant hazard under anticipated conditions of normal use.

Symptoms/effects after inhalation : May cause respiratory irritation.

Symptoms/effects after skin contact : May cause slight irritation.

Symptoms/effects after eye contact : May cause slight irritation. stinging, redness, itching, tears.

Symptoms/effects after ingestion : Ingestion may cause nausea, vomiting and diarrhea. Risk of lung edema.

SECTION 12 Ecological information

12.1. Ecotoxicity

Hazardous to the aquatic environment, short-term (acute) : No data available

Hazardous to the aquatic environment, long-term (chronic) : No data available

12.2. Persistence and degradability

REMOVE Adhesive Remover (bottled)

Persistence and degradability	Not established.
-------------------------------	------------------

C9-12-ISOALKANES (90622-57-4)

Persistence and degradability	Rapidly degradable
-------------------------------	--------------------

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

EUCALYPTUS CITRODORA OIL (8000-48-4)

Persistence and degradability	Rapidly degradable
-------------------------------	--------------------

12.3. Bioaccumulative potential

REMOVE Adhesive Remover (bottled)

Bioaccumulative potential	Not established.
---------------------------	------------------

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

Ozone	: No data available
Fluorinated greenhouse gases	: No
Other information	: Avoid release to the environment.

SECTION 13 Disposal considerations

Product/Packaging disposal recommendations	: Dispose in a safe manner in accordance with local/national regulations.
Additional information	: Handle empty containers with care because residual vapors are flammable.
Ecological waste information	: Avoid release to the environment.

SECTION 14 Transport information

In accordance with DOT / IMDG / IATA

14.1. UN number

UN-No. (DOT)	: UN1993
UN-No. (IMDG)	: 1993
UN-No. (IATA)	: 1993

14.2. UN Proper Shipping Name

Proper Shipping Name (DOT)	: Flammable liquids, n.o.s.(C9-12-ISOALKANES)
Proper Shipping Name (IMDG)	: FLAMMABLE LIQUID, N.O.S. (C9-12-ISOALKANES)
Proper Shipping Name (IATA)	: Flammable liquid, n.o.s.(C9-12-ISOALKANES)

14.3. Transport hazard class(es)

DOT

Transport hazard class(es) (DOT)	: 3
Hazard labels (DOT)	: 3



IMDG

Transport hazard class(es) (IMDG)	: 3
-----------------------------------	-----

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

Hazard labels (IMDG) : 3



IATA
Transport hazard class(es) (IATA) : 3
Hazard labels (IATA) : 3



14.4. Packing group

Packing group (DOT) : III
Packing group (IMDG) : III
Packing group (IATA) : III

14.5. Environmental hazards

Other information : No supplementary information available.

14.6. Transport in bulk

Not applicable

14.7. Special precautions for user

DOT
UN-No. (DOT) : UN1993
DOT Special Provisions (49 CFR 172.102) : B1 - If the material has a flash point at or above 38 C (100 F) and below 93 C (200 F), then the bulk packaging requirements of 173.241 of this subchapter are applicable. If the material has a flash point of less than 38 C (100 F), then the bulk packaging requirements of 173.242 of this subchapter are applicable.
B52 - Notwithstanding the provisions of 173.24b of this subchapter, non-reclosing pressure relief devices are authorized on DOT 57 portable tanks.
IB3 - Authorized IBCs: Metal (31A, 31B and 31N); Rigid plastics (31H1 and 31H2); Composite (31HZ1 and 31HA2, 31HB2, 31HN2, 31HD2 and 31HH2). Additional Requirement: Only liquids with a vapor pressure less than or equal to 110 kPa at 50 C (1.1 bar at 122 F), or 130 kPa at 55 C (1.3 bar at 131 F) are authorized, except for UN2672 (also see Special Provision IP8 in Table 2 for UN2672).
T4 - 2.65 178.274(d)(2) Normal..... 178.275(d)(3)
TP1 - The maximum degree of filling must not exceed the degree of filling determined by the following: Degree of filling = $97 / 1 + a (tr - tf)$ Where: tr is the maximum mean bulk temperature during transport, and tf is the temperature in degrees celsius of the liquid during filling.
TP29 - A portable tank having a minimum test pressure of 1.5 bar (150.0 kPa) may be used provided the calculated test pressure is 1.5 bar or less based on the MAWP of the hazardous materials, as defined in 178.275 of this subchapter, where the test pressure is 1.5 times the MAWP.
DOT Packaging Exceptions (49 CFR 173.xxx) : 150
DOT Packaging Non Bulk (49 CFR 173.xxx) : 203
DOT Packaging Bulk (49 CFR 173.xxx) : 242

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

DOT Quantity Limitations Passenger aircraft/rail (49 CFR 173.27)	: 60 L
DOT Quantity Limitations Cargo aircraft only (49 CFR 175.75)	: 220 L
DOT Vessel Stowage Location	: A - The material may be stowed “on deck” or “under deck” on a cargo vessel and on a passenger vessel.

IMDG	
Special provision (IMDG)	: 223, 274, 955
Limited quantities (IMDG)	: 5 L
Excepted quantities (IMDG)	: E1
Packing instructions (IMDG)	: LP01, P001
IBC packing instructions (IMDG)	: IBC03
Tank instructions (IMDG)	: T4
Tank special provisions (IMDG)	: TP1, TP29
EmS-No. (Fire)	: F-E - FIRE SCHEDULE Echo - NON-WATER-REACTIVE FLAMMABLE LIQUIDS
EmS-No. (Spillage)	: S-E - SPILLAGE SCHEDULE Echo - FLAMMABLE LIQUIDS, FLOATING ON WATER
Stowage category (IMDG)	: A

IATA	
Special provision (IATA)	: A3
PCA Excepted quantities (IATA)	: E1
PCA Limited quantities (IATA)	: Y344
PCA limited quantity max net quantity (IATA)	: 10L
PCA packing instructions (IATA)	: 355
PCA max net quantity (IATA)	: 60L
CAO packing instructions (IATA)	: 366
CAO max net quantity (IATA)	: 220L
ERG code (IATA)	: 3L

SECTION 15 Regulatory information

15.1. Federal regulations

All components of this product are present and listed as Active on the United States Environmental Protection Agency Toxic Substances Control Act (TSCA) inventory, except for:		
C9-12-ISOALKANES	CAS-No. 90622-57-4	≥ 20 – < 50%
Rosemary oil	CAS-No. 84604-14-8	< 0.1%
l-.beta.-Bisabolene	CAS-No. 495-61-4	< 0.1%
ALOE BARBADENSIS LEAF JUICE	CAS-No. 94349-62-9	< 0.1%

Chemical(s) subject to the reporting requirements of Section 313 or Title III of the Superfund Amendments and Reauthorization Act (SARA) of 1986 and 40 CFR Part 372.		
Eugenyl methyl ether	CAS-No. 93-15-2	< 0.1%

n-butyl acetate (123-86-4)	
CERCLA RQ	5000 lb

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

15.2. International regulations

CANADA

C9-12-ISOALKANES (90622-57-4)

Listed on the Canadian DSL (Domestic Substances List)

Dipropylene glycol monomethyl ether (34590-94-8)

Listed on the Canadian DSL (Domestic Substances List)

benzaldehyde (100-52-7)

Listed on the Canadian DSL (Domestic Substances List)

(R)-p-mentha-1,8-diene, d-limonene (5989-27-5)

Listed on the Canadian DSL (Domestic Substances List)

gamma-Undecalactone (104-67-6)

Listed on the Canadian DSL (Domestic Substances List)

Ethyl butyrate (105-54-4)

Listed on the Canadian DSL (Domestic Substances List)

ETHYL METHYLPHENYLGLYCIDATE (77-83-8)

Listed on the Canadian DSL (Domestic Substances List)

ALLYL CAPROATE (123-68-2)

Listed on the Canadian DSL (Domestic Substances List)

1-(2,6,6-TRIMETHYLCYCLOHEX-2-EN-1-YL)PENT-1-EN-3-ONE (1335-46-2)

Listed on the Canadian DSL (Domestic Substances List)

METHOXYPHENYLBUTANONE (104-20-1)

Listed on the Canadian DSL (Domestic Substances List)

CITRAL (5392-40-5)

Listed on the Canadian DSL (Domestic Substances List)

Coumarin (91-64-5)

Listed on the Canadian DSL (Domestic Substances List)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

citronellyl nitrile (51566-62-2)

Listed on the Canadian DSL (Domestic Substances List)

3-P-CUMENYL-2-METHYLPROPIONALDEHYDE (103-95-7)

Listed on the Canadian DSL (Domestic Substances List)

Dimethylcyclohex-3-ene- 1-carbaldehyde (68737-61-1)

Listed on the Canadian DSL (Domestic Substances List)

TRANS-2-HEXENAL (6728-26-3)

Listed on the Canadian DSL (Domestic Substances List)

Linalool (78-70-6)

Listed on the Canadian DSL (Domestic Substances List)

DIMETHYLHEPTENAL (106-72-9)

Listed on the Canadian DSL (Domestic Substances List)

geranyl acetate (105-87-3)

Listed on the Canadian DSL (Domestic Substances List)

ALPHA PINENE NATURAL (80-56-8)

Listed on the Canadian DSL (Domestic Substances List)

beta-Pinene (127-91-3)

Listed on the Canadian DSL (Domestic Substances List)

dl-Limonene (racemic) (138-86-3)

Listed on the Canadian DSL (Domestic Substances List)

BETA-CARYOPHYLLENE (87-44-5)

Listed on the Canadian DSL (Domestic Substances List)

Myrcene (123-35-3)

Listed on the Canadian DSL (Domestic Substances List)

CAMPHENE (79-92-5)

Listed on the Canadian DSL (Domestic Substances List)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

Longifolene (475-20-7)

Listed on the Canadian DSL (Domestic Substances List)

Terpinolene (586-62-9)

Listed on the Canadian DSL (Domestic Substances List)

METHYL SALICYLATE (119-36-8)

Listed on the Canadian DSL (Domestic Substances List)

DECANAL (112-31-2)

Listed on the Canadian DSL (Domestic Substances List)

3,7-Диметилокта-1,6-диен-3-ол (78-70-6)

Listed on the Canadian DSL (Domestic Substances List)

Alpha-Terpineol (98-55-5)

Listed on the Canadian DSL (Domestic Substances List)

EUCALYPTOL (470-82-6)

Listed on the Canadian DSL (Domestic Substances List)

Kamfer; Camphor, synthetic; Камфор; Kamferi; Bornan-2-on; Bornan-2-one; Kamfóra; Kamparas; Kamfer (syntetisk); Kamfora syntetyczna bornan-2-on; Cànfora (sintética); Camfor; kafor; 1,7,7-Триметилбицикло [2,2,1]гептан-2-он; Alcanfor sintético; Campher (s. Kamfer); Alcator sintético; Camphre; Camphre (synthétique) # Kamfer (synthetisch); Kamfor, tiruan # Camphor, synthetic; Gáfor (bornán-2-ón); Campher, syntetisk (76-22-2)

Listed on the Canadian DSL (Domestic Substances List)

(S)-p-mentha-1,8-diene, l-limonene (5989-54-8)

Listed on the Canadian DSL (Domestic Substances List)

L-BORNEOL (464-45-9)

Listed on the Canadian DSL (Domestic Substances List)

ALPHA-TERPINEOL (98-55-5)

Listed on the Canadian DSL (Domestic Substances List)

p-Cymene (99-87-6)

Listed on the Canadian DSL (Domestic Substances List)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

laevo-Bornyl acetate (5655-61-8)

Listed on the Canadian DSL (Domestic Substances List)

Linalool (78-70-6)

Listed on the Canadian DSL (Domestic Substances List)

3-CARENE (13466-78-9)

Listed on the Canadian DSL (Domestic Substances List)

l-.beta.-Bisabolene (495-61-4)

Not listed on the Canadian DSL (Domestic Substances List)/NDSL (Non-Domestic Substances List)

Oct-1-ene-3-ol (3391-86-4)

Listed on the Canadian DSL (Domestic Substances List)

cis-beta-Ocimene (3338-55-4)

Listed on the Canadian DSL (Domestic Substances List)

Eugenyl methyl ether (93-15-2)

Listed on the Canadian DSL (Domestic Substances List)

Cinnamaldehyde (104-55-2)

Listed on the Canadian DSL (Domestic Substances List)

Citronellal (106-23-0)

Listed on the Canadian DSL (Domestic Substances List)

3,7-dimethyl-6-octen-1-ol (citronellol) (106-22-9)

Listed on the Canadian DSL (Domestic Substances List)

3,7-Dimethyl-6-octen-1-ol acetate ; Citronellyl acetate (150-84-5)

Listed on the Canadian DSL (Domestic Substances List)

Isopulegol (89-79-2)

Listed on the Canadian DSL (Domestic Substances List)

Geranyl acetate (105-87-3)

Listed on the Canadian DSL (Domestic Substances List)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

CINNAMON LEAF OIL (8015-91-6)
Listed on the Canadian DSL (Domestic Substances List)
n-butyl acetate (123-86-4)
Listed on the Canadian DSL (Domestic Substances List)
isoamyl salicylate (87-20-7)
Listed on the Canadian DSL (Domestic Substances List)
COUMARIN (91-64-5)
Listed on the Canadian DSL (Domestic Substances List)
CEDARWOOD OIL (8000-27-9)
Listed on the Canadian DSL (Domestic Substances List)
LINALYL ACETATE (115-95-7)
Listed on the Canadian DSL (Domestic Substances List)
Diphenylether;Phenyl ether (vapour);Дифенилетер (Дифенилоксид);苯基醚(二苯醚) # Phenyl ether;Difenylether;Difenüüleeter;Fenyylieetteri;Difenil-eter;Diphenyl ether (vapour);Dífenýleter (dífenýloxíd, fenýleter);Difenyleter (Fenyleter);Eter difenylowy;Éter fenílico (vapor);Fenileter (vapor);difenileter – para;Diphenyl ether;Phenyl ether, vapor;Оксидибензол;Feniléter (Éter fenílico);Feniléter;Biphenylether (s. Diphenylether);Diphenylether (Dampf);Eter fenílico, vapores;Oxyde de biphényle;Oxyde de diphényle (vapeur) # Difenyloxide (damp);Fenil eter, wap # Phenyl ether, vapour;Difenyléter;Diphenylether (Diphenyloxid; Phenylether) (101-84-8)
Listed on the Canadian DSL (Domestic Substances List)
Linalool (78-70-6)
Listed on the Canadian DSL (Domestic Substances List)
HELIOTROPINE (120-57-0)
Listed on the Canadian DSL (Domestic Substances List)
TRANS-ANETHOLE (4180-23-8)
Listed on the Canadian DSL (Domestic Substances List)
GLYCINE SOJA OIL (8001-22-7)
Listed on the Canadian DSL (Domestic Substances List)
ALOE BARBADENSIS LEAF JUICE (94349-62-9)
Not listed on the Canadian DSL (Domestic Substances List)/NDSL (Non-Domestic Substances List)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

Benzyl alcohol (100-51-6)

Listed on the Canadian DSL (Domestic Substances List)

EU-Regulations

No additional information available

National regulations

C9-12-ISOALKANES (90622-57-4)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Dipropylene glycol monomethyl ether (34590-94-8)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

benzaldehyde (100-52-7)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

(R)-p-mentha-1,8-diene, d-limonene (5989-27-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

gamma-Undecalactone (104-67-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Ethyl butyrate (105-54-4)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

1-(2,6,6-TRIMETHYLCYCLOHEX-2-EN-1-YL)PENT-1-EN-3-ONE (1335-46-2)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

METHOXYPHENYLBUTANONE (104-20-1)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

CITRAL (5392-40-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Coumarin (91-64-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

citronellyl nitrile (51566-62-2)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

TRANS-2-HEXENAL (6728-26-3)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Linalool (78-70-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

DIMETHYLHEPTENAL (106-72-9)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

geranyl acetate (105-87-3)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

ALPHA PINENE NATURAL (80-56-8)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

beta-Pinene (127-91-3)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

dl-Limonene (racemic) (138-86-3)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Myrcene (123-35-3)

Listed on IARC (International Agency for Research on Cancer)
Listed on INSQ (Mexican National Inventory of Chemical Substances)

CAMPHERE (79-92-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Terpinolene (586-62-9)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

METHYL SALICYLATE (119-36-8)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

DECANAL (112-31-2)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

3,7-Диметилокта-1,6-диен-3-ол (78-70-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

Alpha-Terpineol (98-55-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

EUCALYPTOL (470-82-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Kampfer; Camphor, synthetic; Камфор; Kamferi; Bornan-2-on; Bornan-2-one; Kamfóra; Kamparas; Kamfer (syntetisk); Kamfora syntetyczna bornan-2-on; Cànfora (sintética); Camfor; kafora; 1,7,7-Триметилбицикло [2,2,1]гептан-2-он; Alcanfor sintético; Campher (s. Kampfer); Alcator sintético; Camphre; Camphre (synthétique) # Kamfer (synthetisch); Kamfor, tiruan # Camphor, synthetic; Gáfor (bornán-2-ón); Campher, syntetisk (76-22-2)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

ALPHA-TERPINEOL (98-55-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

p-Cymene (99-87-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Linalool (78-70-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Oct-1-ene-3-ol (3391-86-4)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

cis-beta-Ocimene (3338-55-4)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Eugenyl methyl ether (93-15-2)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)
Listed on INSQ (Mexican National Inventory of Chemical Substances)

Cinnamaldehyde (104-55-2)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Citronellal (106-23-0)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

3,7-dimethyl-6-octen-1-ol (citronellol) (106-22-9)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

3,7-Dimethyl-6-octen-1-ol acetate ; Citronellyl acetate (150-84-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Geranyl acetate (105-87-3)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

n-butyl acetate (123-86-4)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

COUMARIN (91-64-5)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

LINALYL ACETATE (115-95-7)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Diphenylether;Phenyl ether (vapour);ДифенилETER (Дифенилоксид);苯基醚(二苯醚) # Phenyl ether;Difenylether;Difenüüleeter;Fenyylieetteri;Difenil-eter;Diphenyl ether (vapour);Dífenýleter (dífenýloxíd, fenýleter);Difenyleter (Fenyleter);Eter difenylovy;Éter fenílico (vapor);Fenileter (vapor);difenileter – para;Diphenyl ether;Phenyl ether, vapor;Оксидибензол;Feniléter (Éter fenílico);Feniléter;Biphenylether (s. Diphenylether);Diphenylether (Dampf);Eter fenílico, vapores;Oxyde de biphényle;Oxyde de diphényle (vapeur) # Difenyloxide (damp);Fenil eter, wap # Phenyl ether, vapour;Difenyléter;Diphenylether (Diphenyloxid; Phenylether) (101-84-8)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

Linalool (78-70-6)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

HELIOTROPINE (120-57-0)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

TRANS-ANETHOLE (4180-23-8)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

GLYCINE SOJA OIL (8001-22-7)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

ALOE BARBADENSIS LEAF JUICE (94349-62-9)

Listed on INSQ (Mexican National Inventory of Chemical Substances)

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

Benzyl alcohol (100-51-6)
Listed on INSQ (Mexican National Inventory of Chemical Substances)

15.3. State regulations

 **WARNING:** This product can expose you to beta-Myrcene, which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

Component	State or local regulations
Dipropylene glycol monomethyl ether(34590-94-8)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List
benzaldehyde(100-52-7)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List
Ethyl butyrate(105-54-4)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List
ALPHA PINENE NATURAL(80-56-8)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List
dl-Limonene (racemic)(138-86-3)	U.S. - New Jersey - Right to Know Hazardous Substance List
Terpinolene(586-62-9)	U.S. - New Jersey - Right to Know Hazardous Substance List
METHYL SALICYLATE(119-36-8)	U.S. - Pennsylvania - RTK (Right to Know) List
Kampfer;Camphor, synthetic;Камфор;Kamferi;Bornan-2-on;Bornan-2-one;Kamfóra;Kamparas;Kamfer (syntetisk);Kamfora syntetyczna bornan-2-on;Cânfora (sintética);Camfor;kafra;1,7,7-Триметилбицикло [2,2,1]гептан-2-он;Alcanfor sintético;Campher (s. Kampfer);Alcantor sintético;Camphre;Camphre (synthétique) # Kamfer (synthetisch);Kamfor, tiruan # Camphor, synthetic;Gáfor (bornán-2-ón);Campher, syntetisk(76-22-2)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List
p-Cymene(99-87-6)	U.S. - Massachusetts - Right To Know List; U.S. - Pennsylvania - RTK (Right to Know) List
Eugenyl methyl ether(93-15-2)	U.S. - New Jersey - Right to Know Hazardous Substance List
n-butyl acetate(123-86-4)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List

REMOVE Adhesive Remover (bottled)

Safety Data Sheet

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)

Component	State or local regulations
Diphenylether;Phenyl ether (vapour);Дифенилетер (Дифенилоксид);苯基醚(二苯醚) # Phenyl ether;Difenylether;Difenüüleeter;Fenyylieetteri;Difenil-eter;Diphenyl ether (vapour);Dífenýleter (dífenýloxíd, fenýleter);Difenyleter (Fenyleter);Eter difenylowy;Éter fenílico (vapor);Fenileter (vapori);difenileter – para;Diphenyl ether;Phenyl ether, vapor;Оксидибензол;Feniléter (Éter fenílico);Feniléter;Biphenylether (s. Diphenylether);Diphenylether (Dampf);Eter fenílico, vapores;Oxyde de biphényle;Oxyde de diphényle (vapeur) # Difenyloxide (damp);Fenil eter, wap # Phenyl ether, vapour;Dífenýleter;Diphenylether (Diphenyloxid; Phenylether)(101-84-8)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List
GLYCINE SOJA OIL(8001-22-7)	U.S. - Pennsylvania - RTK (Right to Know) List
Benzyl alcohol(100-51-6)	U.S. - Massachusetts - Right To Know List; U.S. - Pennsylvania - RTK (Right to Know) List

SECTION 16 Other information

according to 29 CFR § 1910.1200, Hazard Communication Standard (HCS 2024)	
Issue date	: 8/7/2026
Data sources	: REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labeling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.
Other information	: None.

Full text of hazard classes and H-statements	
H226	Flammable liquid and vapor
H304	May be fatal if swallowed and enters airways
H315	Causes skin irritation
H317	May cause an allergic skin reaction
H319	Causes serious eye irritation
H411	Toxic to aquatic life with long lasting effects

Safety Data Sheet (SDS), USA

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.