


 **smith&nephew**

TWINFIX[®] Ultra

Preloaded Suture Anchors
with Needles



Needle management with a twist

TWINFIX[◊] Ultra Preloaded Suture Anchors with Needles

Easy-to-use handles for speed and efficiency. TWINFIX Ultra Preloaded Suture Anchors with Needles combine best-in-class anchor pullout strength¹ with controlled suture and needle management for your mini-open rotator cuff repairs.

Trusted anchor material choices include absorbable (HA), non-absorbable (PK) and metallic (Ti) materials.

Increased visibility; preloaded blue suture brings reliable suture fixation with low knot profile.

Secure fixation is achieved with a transitional thread design. Distal “cutting” threads provide easy insertion and proximal “locking” threads provide strong pullout strength, even in poor bone quality.¹

Preloaded and ready-to-use needles in handle.

Ergonomically designed handles with twist-and-release needle delivery system.


After placing anchor in the bone, please follow the steps below:



Step 1 - Ready



Step 2 - Twist



Step 3 - Full Turn

Ordering Information

TWINFIX Ultra with Needles available in 4.5 mm, 5.5 mm and 6.5 mm sizes

Reference #	Description
72202611	TWINFIX Ultra PK 4.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202615	TWINFIX Ultra PK 5.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202619	TWINFIX Ultra PK 6.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202612	TWINFIX Ultra PLLA/HA 4.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202616	TWINFIX Ultra PLLA/HA 5.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202620	TWINFIX Ultra PLLA/HA 6.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202613	TWINFIX Ultra Ti 4.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202618	TWINFIX Ultra Ti 5.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles
72202622	TWINFIX Ultra Ti 6.5 mm Suture Anchor with two #2 ULTRABRAID [◊] Sutures (blue, blue-cobraid) with Needles

Endoscopy
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810 USA

www.smith-nephew.com
T +1 978 749 1000
US Customer Service: +1 800 343 5717
International Customer Service: +1 978 749 1140

¹Data on file at Smith & Nephew 2009, ITR-4091, ITR-4039, ITR-3952.

[◊]Trademark of Smith & Nephew. Certain marks Reg. US Patent and Trademark Office. All Trademarks acknowledged. ©2011 Smith & Nephew. All rights reserved. Printed in USA. 06/11 10600758 Rev. A