

# + Elevate efficiency

Hip Repair Solutions


Smith+Nephew



CAP-LIFT<sup>◇</sup>  
Cannula



## Designed to transform with a unique open style



### Increases procedural efficiency<sup>1</sup>

- Enables hands-free, **customizable capsular elevation**<sup>1</sup>
- Increases visibility and motion<sup>\*1</sup>
- Eliminates switching sticks and sleds<sup>1</sup>
- **Reduces procedural steps** by simplifying portal changes and instrument exchanges<sup>\*\*1</sup>
- **Design to prevent iatrogenic cartilage damage** with capsule lock feature

<sup>\*1</sup>When compared to standard cannulas

<sup>\*\*1</sup>When compared to intermediate instrument switching steps and devices

## Ordering information

CAP-LIFT <sup>®</sup> Cannula	
Reference	Description
72205889	CAP-LIFT Cannula, 70mm
72205912	CAP-LIFT Cannula, 80mm
72205890	CAP-LIFT Cannula, 90mm
72205913	CAP-LIFT Cannula, 100mm
72205891	CAP-LIFT Cannula, 110mm

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

## Learn more at [smith-nephew.com](https://smith-nephew.com)

Smith & Nephew, Inc.  
150 Minuteman Road  
Andover, MA 01810  
USA

[www.smith-nephew.com](https://www.smith-nephew.com)  
T +978 749 1000  
US Customer Service:  
+1 800 343 5717

<sup>®</sup>Trademark of Smith+Nephew.  
All trademarks acknowledged.  
©2025 Smith+Nephew. All rights  
reserved. 46801 V3 11/25

### References

1. Smith+Nephew Internal Report 10151735 Rev A