

# + The future in your hands

Backed by 20 years of experience,  
legacy and innovation, get ready to  
experience the future in your hands

**Smith+Nephew**

FLOW 90<sup>◇</sup>  
COBLATION<sup>®</sup> Wand





FLOW 90° has been shown to have a significantly faster ablation rate than other frequently used 90-degree arthroscopy wands<sup>4\*</sup>

## A new innovative tip design...

- Triangular tip is designed to enable tissue manipulation
- Fits down a 5.0mm cannula
- Electrode engineered to tackle complex pathology

## ...combined with proven technology...

- FLOW~IQ° technology changes the size and intensity of the plasma field based on desired tissue effects to address various shoulder pathologies
- AMBIENT° technology provides accurate, real-time temperature monitoring of the intra-articular fluid<sup>1,2</sup>
- SCOPE-SENSING° technology provides a higher level of protection for expensive surgical equipment<sup>3\*\*</sup>
- FLOW 90's integrated suction tubing eliminates the cost associated with a separate suction tubing line

## Ordering information

FLOW 90° COBLATION° wand	
Reference #	Description
72290146	WEREWOLF+ Generator
72290038	FLOW 90 wand

\*As demonstrated *in vitro* when tested at default settings;  $p < 0.05$

\*\*compared to QUANTUM° 2, demonstrated *in vitro*

Learn more at [werewolfcoblation.com](http://werewolfcoblation.com)

Distributed by:  
Sports Medicine  
Smith & Nephew, Inc.  
150 Minuteman Road  
Andover, MA 01810

Manufactured by:  
ArthroCare Corporation  
7000 West William Cannon Drive  
Austin, TX 78735, USA

[www.smith-nephew.com](http://www.smith-nephew.com)  
T +978 749 1000  
US Customer Service:  
+1 800 343 5717

°Trademark of Smith & Nephew.  
©2022 Smith & Nephew. All  
rights reserved. All trademarks  
acknowledged. 17186 V3 02/22

### References

1. ArthroCare Corporation 2015. Report, Design Verification, FLOW 50 at Maximum 1X Sterility. P/N 65200-01 Rev.A. 2. Smith+Nephew 2019. Report, Design Verification, 1X Ablation Testing & Every Tissue Every Mode (ETEM), FLOW 90. P/N 88079-01 Rev.C. 3. Smith+Nephew 2019. Report, SCOPE-SENSING Efficacy Study, Quantum 2 and Werewolf. P/N 106506-01 Rev.A. 4. Smith+Nephew 2020. Report, Performance Characterization, 90 Degree Sports Medicine Wands Manufactured by Smith & Nephew, DePuy, Stryker, Arthrex, and Conmed. P/N 112090 Rev.A.