



# Extend the capabilities

of ENDOBUTTON® CL ULTRA Fixed Loop and  
ULTRABUTTON® Adjustable Loop Fixation Devices

 **smith&nephew**  
**XTENDOBUTTON®**  
Family of Fixation Devices

Supporting healthcare professionals

# Reconstruct with confidence


## Enhance your cortical fixation options

**XTENDOBUTTON® Oval Fixation Devices** are designed to nestle onto native anatomy



**Peg design**  
Centers button in the tunnel


**XTENDOBUTTON® Round Fixation Devices** are designed to provide extended surface area



**Suture management holes**  
Option to tie additional sutures

**Extended design**  
Distributes load over greater surface area than the oval

Easy to connect without removing graft



### Ordering information

XTENDOBUTTON Round Fixation Devices		
Reference #	Tunnel diameter	Description
71935601	4.0mm-7.0mm	XTENDOBUTTON Round 15mm ø (small)
71935602	7.5mm-9.5mm	XTENDOBUTTON Round 17mm ø, 7.5mm peg (medium)
71935603	10.0mm-13.0mm	XTENDOBUTTON Round 21mm ø, 10mm peg (large)
XTENDOBUTTON Oval Fixation Devices		
Reference #	Tunnel diameter	Description
71935604	4.0mm-7.0mm	XTENDOBUTTON Oval 9mm x 16mm (small)
71935605	7.5mm-9.5mm	XTENDOBUTTON Oval 11mm x 18mm, 7.5mm peg (medium)
71935606	10.0mm-13.0mm	XTENDOBUTTON Oval 15mm x 21mm, 10mm peg (large)

XTENDOBUTTON® Fixation Device		
Reference #	Tunnel diameter	Description
72200134	6.0-10.0mm	XTENDOBUTTON Fixation Device

ULTRABUTTON® Adjustable Fixation Device	
Reference #	Description
72290003	ULTRABUTTON Adjustable

ENDOBUTTON® CL ULTRA Fixation Devices	
Reference #	Description
772203331	10mm ENDOBUTTON CL ULTRA
72200146-72200155	15mm-60mm ENDOBUTTON CL ULTRA (5mm increments)

To order these products, call **+1 800 343 5717** in the U.S. or contact an authorized Smith & Nephew representative.

**Sports Medicine**  
Smith & Nephew, Inc.  
150 Minuteman Road  
Andover, MA 01810

[www.smith-nephew.com](http://www.smith-nephew.com)  
T +978 749 1000  
US Customer Service:  
+1 800 343 5717

®Trademark of Smith & Nephew.  
©2019 Smith & Nephew.  
All rights reserved. Printed  
in USA. 14416 V2 06/19

The ULTRABUTTON Adjustable Fixation Device is manufactured by ArthroCare Corporation, 7000 West William Cannon Drive, Austin, TX 78735, USA. All other products listed here are manufactured by Smith & Nephew.