# The future in your hands

## Introducing the new FLOW 90 Wand


Backed by 20 years of experience, legacy and innovation, get ready to experience the future in your hands.

## Smith-Nephew

WEREWOLF FLOW 90 Wand with FLOW -IQ Technology







#### ...combined with proven technology...

- FLOW~IQ<sup>o</sup> technology changes the size and intensity of the plasma field based on desired tissue effects to address various shoulder pathologies
- AMBIENT<sup>o</sup> technology provides accurate, real-time temperature monitoring of the intra-articular fluid<sup>1,2</sup>
- SCOPE-SENSING<sup>o</sup> technology provides a higher level of protection for expensive surgical equipment<sup>3</sup>\*
- FLOW 90's integrated suction tubing eliminates the cost associated with a separate suction tubing line

### Ordering information

WEREWOLF° COBLATION° FLOW 90° WAND	
Reference #	Description
72290105	WEREWOLF Generator
72290038	FLOW 90 wand

<sup>\*</sup>compared to QUANTUM<sup> 2</sup>, demonstrated in vitro

Smith & Nephew Pty Ltd Smith & Nephew Ltd www.smith-nephew.com/australia Australia New Zealand

T+61 2 9857 3999 T +64 9 820 2840 F+61 2 9857 3900 F+64 9 820 2841 www.smith-nephew.com/new-zealand

indications for use, contraindications, precautions <sup>o</sup>Trademark of Smith & Nephew. All trademarks and warnings, please consult the product's acknowledged. 17186-anz V2 SN14860 REV0 03/20 applicable Instructions for Use (IFU) prior to use.

For detailed product information, including