



Smith+Nephew

ACUTE QUICK
CONNECT STRUTS
Circular Fixator




- + Poly-axial Struts for rapid circular frame assembly
- + A typical Pilon frame consists of 3 Pins in the tibia and 3 Wires in the foot
- + Four Struts are required for a weight-bearing construct





Cat. No.	Description	Collapsed	Extended	Range
71935866	ACUTE QC STRUT XShort	125mm	150mm	25mm
71935867	ACUTE QC STRUT Short	140mm	180mm	40mm
71935868	ACUTE QC STRUT Medium	175mm	255mm	80mm
71935869	ACUTE QC STRUT Long	230mm	350mm	120mm
71935870	Speed Nut	sold in pack of 8		
71935871	ACUTE QC STRUT Band	sold in pack of 4		



[Watch a surgical demonstration](#) of Pilon fracture stabilization using ACUTE STRUTS by Prof Nando Ferreira

Prof Nando Ferreira
Stellenbosch University, Cape Town

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, Tennessee 38116
USA

www.smith-nephew.com

°Trademark of Smith+Nephew
All Trademarks acknowledged
©2024 Smith+Nephew
45563 V1 10/24