


# 絆創膏に ワンランク上のやさしさを

浅い傷



皮膚裂傷  
(スキん-テア)



保護



smith&nephew

ハイドロ  
ジェントルエイド

Supporting healthcare professionals  
for over 150 years

## 傷を優しく守る

クッション性に優れたふんわりパッドが肌によくなじみ、傷を守ります。

## 最長7日間\*

ポリウレタンフォームが滲出液を吸収し、最長で7日間継続して使用できます。[3]

イメージ図

過剰な  
滲出液を素早く蒸散


細菌の侵入を  
防ぎます

滲出液を  
吸収・保持

\*創を観察し、滲出液の量に応じて適宜交換をしてください。

## 防水仕様


水仕事やシャワー浴でも剥がれにくい防水フィルムを採用しました。[1][2]



## 肌に優しい

シリコン粘着剤がソフトに肌に密着し、剥がす時の痛みを軽減します。[3][4]  
敏感肌にも優しい低アレルギー粘着剤です。[5][6][7]

イメージ図



従来の  
アクリル系粘着剤

肌に優しい  
一般的な粘着剤

シリコン  
粘着剤

粘着力が落ちにくいいため、ずれた時などに貼り直しができます。

(感じ方には個人差があります)

## 使用例

浅い傷・保護



小児の  
すり傷/切傷



# ハイドロ ジェントルエイド


一般医療機器 救急絆創膏  
届出番号:  
13B1X10222WC0016

製品番号	サイズ	パッド部サイズ	枚/箱	標準価格									
66800833	5×5 cm	3×3 cm	10	¥1,590									
66800834	7.5×7.5 cm	5×5 cm	10	¥1,950									
66800835	10×10 cm	7.5×7.5 cm	10	¥2,400									
66800840	15×15 cm	12.5×12.5 cm	10	¥5,390									
66800836	5.5×12 cm	10	¥2,330	66800837	8×15 cm	5.5×12.5 cm	10	¥3,570	66801012	10×20 cm	7.2×17.2 cm	10	¥4,340
66800837	8×15 cm	5.5×12.5 cm	10	¥3,570									
66801012	10×20 cm	7.2×17.2 cm	10	¥4,340									

**使用目的又は効果** 傷及び周囲皮膚の被覆及び保護

**禁忌・禁止\*** 再使用禁止

\*詳細は添付文書をご参照ください。



5×5 cm  
(66800833)

7.5×7.5 cm  
(66800834)


10×10 cm  
(66800835)

15×15 cm  
(66800840)

5.5×12 cm  
(66800836)

8×15 cm  
(66800837)

10×20 cm  
(66801012)



## Reference

- [1] L. Tompkins. ALLEVYN Gentle Border Lite Dressing Physical Properties. Smith & Nephew Medical Ltd, Hull DS/10/100/R1, Aug 2010.
- [2] H. Dharma. An open, prospective, randomized, within volunteer, comparisons of ALLEVYN Gentle Border and Mepilex Border. Smith & Nephew Medical Ltd, Hull OR-DOF/011, 24 Aug 2010.
- [3] Hurd, T et al. A multi-centre in-market evaluation of ALLEVYN Gentle Border. Wounds UK 2009;5(3):32-44.
- [4] K. Vowden et al. An open prospective randomized, multi-centre clinical evaluation of a hydrocellular silicone foam dressing in the management of exuding chronic and acute wounds. Presented at the EWMA, Brussels 2011.
- [5] Stephen Haynes, J et al. Clinical evaluation of ALLEVYN Gentle Border Lite within one health trust. BJN;2011;20:S36,S38-42.
- [6] Denyer JE. Wound management for children with epidermolysis bullosa. Dermatol Clin. 2010 Apr;28(2):257-64.
- [7] Walton G. Product Safety Statement 196"ALLEVYN Gentle Border Lite", Smith & Nephew Ltd, Hull 2016

## 製造販売元

スミス・アンド・ネフュー株式会社 ウンドマネジメント事業部

〒105-0011 東京都港区芝公園二丁目4番1号 TEL.03-5403-8830

札幌営業所 TEL.011-736-9596 FAX.011-736-9755

仙台営業所 TEL.022-276-6726 FAX.022-276-6652

東京営業所 TEL.03-5403-8915 FAX.03-5403-8916

名古屋営業所 TEL.052-221-1407 FAX.052-221-1420

大阪営業所 TEL.06-6399-3303 FAX.06-6399-3993

広島営業所 TEL.082-502-7322 FAX.082-502-7325

福岡営業所 TEL.092-452-0142 FAX.092-452-0146

<http://www.smith-nephew.com/japan/>

© 2019 Smith & Nephew KK

## 販売店

