



Control the **flow** for optimal performance on multiple tissue types

 **smith&nephew**

**COBLATION<sup>◇</sup>**

**FLOW 50<sup>◇</sup>**

Wand with FLOW~IQ<sup>◇</sup> Technology


Supporting healthcare professionals

# Control the **flow** for optimal performance on multiple tissue types

- FLOW~IQ<sup>®</sup> Technology provides automated outflow control that delivers a new standard of ablation performance.
- The FLOW 50<sup>°</sup> Wand adapts to multiple tissue types such as articular cartilage, meniscus, and ligamentous tissue, efficiently resecting targeted tissue while preserving healthy surrounding tissue.



Improved Integrated Finger Switch (IFS) Technology.



Slim electrode design facilitates access to targeted tissue.



Cartilage treatment



Meniscal sculpting



ACL debridement

## Technical specifications

Catalog number	72290037
Controller	WEREWOLF Controller
Shaft Diameter	3.76 mm (0.148")
Working length	135 mm (5.3")
Tip Width	2.95 mm (0.116")
Default set point	Med

## Ordering information

Product Name	
Reference #	Description
72290037	FLOW 50 Wand

ArthroCare Corporation  
7000 West William Cannon Drive  
Austin, Texas 78735  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)  
T +978 749 1000  
US Customer Service:  
+1 800 343 5717

<sup>®</sup>Trademark of Smith & Nephew.  
©2019 Smith & Nephew.  
All rights reserved. Printed in USA.  
19122 V1 P/N 68924 Rev. A 04/19



Supporting healthcare professionals for over 150 years