

# + PROFORE<sup>◇</sup> Multi-layer Compression Bandage System

Delivering up to 7 days sustained,  
graduated pressure<sup>1-4</sup>



## Smith+Nephew

PROFORE<sup>◇</sup>  
Multi-layer Compression  
Bandage System

[www.smith-nephew.com/profore](http://www.smith-nephew.com/profore)

Ask your pharmacist or doctor  
about PROFORE

# The **PROFORE** Multi-layer Compression Bandage System is a range of multi-layered compression bandage systems that provide sustained graduated pressure<sup>1-3</sup> for up to 7 days<sup>1,2,4</sup>

## Features and benefits



### Sustained graduated pressure

The multi-layered bandage technique provides sustained graduated pressure.<sup>1,2,3</sup> The PROFORE Bandage System is designed to provide higher pressure at the ankle with decreasing pressure at the calf<sup>1,2</sup>



### Up to 7 day wear time

The PROFORE Bandage System has been demonstrated to maintain effective levels of compression<sup>1</sup> whilst being able to remain in place for up to 7 days<sup>1,2,4</sup>



### Provides recommended levels of compression

The PROFORE Bandage System has been clinically proven to provide recommended levels of compression of 40mmHg pressure at the ankle decreasing to 17mmHg just below the knee<sup>\*1,2</sup>



### Demonstrated clinical efficacy

The PROFORE Bandage System has been shown to achieve complete wound closure in up to 72% of venous ulcers by 12 weeks and up to 84% by 24 weeks<sup>\*5</sup>



### Patient comfort

The padding bandage layer (PROFORE #1) is soft and conformable to support ease of application.<sup>6</sup> The absorbent properties of the bandage system reduce the risk of skin maceration and irritation caused by excess wound exudate<sup>6</sup>



### Range of options including latex free and reduced compression

The PROFORE Bandage System range includes latex free formulations for those patients with latex allergies<sup>7</sup> and reduced compression versions (LITE) developed specially for the management of mixed aetiology leg ulcers<sup>8,9</sup>

\*Excluding PROFORE LITE.

Advanced Wound Management  
Smith & Nephew Medical Ltd  
101 Hessle Road  
Hull HU3 2BN  
UK  
T +44 (0) 1482 225181  
F +44 (0) 1482 328326

[www.smith-nephew.com](http://www.smith-nephew.com)

©Trademark of Smith+Nephew  
All Trademarks acknowledged  
©March 2022 Smith+Nephew  
AWM-AWD-16364. GMC1489

## References

1. Blair SD, Wright D, Backhouse CM, Riddle E, Mccollum CN. Sustained compression and healing of chronic venous ulcers. *BMJ*. 1988;297(6657):1159–1161. 2. Moffatt C, Stubbings N, Richardson B. The Charing Cross approach to venous ulcers: Christine Moffatt and Nicola Stubbings examine the challenge that venous ulcers pose and one strategy for dealing with them. *Nursing Standard*. 1990;5(12):6–13. 3. Ukat A, Konig M, Vanscheidt W, Münter K. Short-stretch versus multilayer compression for venous leg ulcers: a comparison of healing rates. *Journal of wound care*. 2003;12(4):139–143. 4. Bowering CK. Use of Layered Compression Bandages in Diabetic Patients. *Advances in Wound Care*. 1998;11:129–135. 5. Ukat A, Konig M, Vanscheidt W, Münter K. Short-stretch versus multilayer compression for venous leg ulcers: a comparison of healing rates. *Journal of wound care*. 2003;12(4):139–143. 6. Donald A. 19 Dec 1994. Soffban Natural Questionnaire. Internal report. 7. Smith+Nephew 2018. Regulatory Compliance Document - Raw Materials Profore#1. Internal report. 8. Kwon C, Hegarty M, Oxenham W, et al. An indirect testing approach for characterizing pressure profiles of compression bandages and hosiery. *The Journal of The Textile Institute*. 2018;109(2):256–267. 9. Mosti G, Iabichella ML, Partsch H. Compression therapy in mixed ulcers increases venous output and arterial perfusion. *Journal of vascular surgery*. 2012;55(1):122–128. 10. DS/16/373/R Wound model of ALLEVYN Life Non-Bordered vs Mepilex™ XT dressings - 7 day vertical testing on a curved plate under compression bandaging with horse serum using a flow rate modelling that of a moderately exuding wound. October 2016. 11. Jennings S. (2019) ALLEVYN Gentle 2 - Vertical wound model using horse serum at a flow rate modelling that of a moderately exuding wound, using a compression bandage system. Data on file DS/19/197/R.

## Indications

The PROFORE Multi-layer Compression Bandage System has been specifically developed for the management of venous leg ulcers and associated conditions.

An ideal combination for use with ALLEVYN® LIFE NON BORDERED Foam Dressing and ALLEVYN GENTLE Foam Dressing.<sup>10,11</sup>



## PROFORE System

S+N Code	Description	Carton
66000017	25cm - 30cm	4
66000018	Greater than 30cm	8

## PROFORE System components

S+N Code	Description	Carton
66000415	PROFORE LITE	8
66000333	PROFORE # 1	12
66000328	PROFORE # 2	12
66000334	PROFORE # 3	12
66000336	PROFORE # 4	12
66000335	PROFORE +	12
66000701	WCL	50

## PROFORE Latex Free formulation system

S+N Code	Description	Carton
66000770	18cm - 25cm	8
66000771	PROFORE LITE	8

## PROFORE Latex Free formulation system components

S+N Code	Description	Carton
66000757	PROFORE # 1	12
66000758	PROFORE # 2	12
66000778	PROFORE # 3	12
66000779	PROFORE # 4	12
66000776	PROFORE +	12
66000701	WCL	50

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.