



UNIQUE NEEDS, **COMPLETE CARE.**



ALLEVYN[◇]
COMPLETE CARE
Foam Dressing

Smith+Nephew



ALLEVYN[®] COMPLETE CARE Dressing

Every wound tells a different story

And behind every dressing change, there is a patient living with hope – relying on you to help make a difference.

When estimates suggest that up to **1 in 50 adults**¹ are affected by chronic wounds in their lifetime, and **1 out of 10** hospitalised patients are impacted by HAPIs,² confidence in your choice of foam dressing matters now more than ever.

It's a decision that has the potential to change lives

With over 35 years of foam dressing innovation, we have a deep expertise and a proud heritage in advanced wound management. We understand the challenges you face – and the need to make every moment of care count.

We recognise the need to help you do your job better, with a purpose-built solution that helps deliver performance

Because we know wounds – and you know your patients

Together, we have the opportunity to combine innovation and insight. To help progress wound healing, protect fragile skin, and get patients back to what they love.

Your choice has never been so important

Every wound tells a story...

... let's help patients like Zoe* rewrite the outcome

*Zoe does not represent an actual patient

Because every wound and every patient deserves complete care




Patients like Zoe* can get back to who they love with less concerns

about disruption to her day from leakage, detachment or dressing changes

ALLEVYN[◇] COMPLETE CARE Foam Dressing is a true next generation dressing

Designed for strong performance in both pressure injury prevention and wound management, it's built to meet you daily challenges and help support healing.^{3,4}

Most importantly, it's built to get patients like Zoe **back to living the moments they love.**



By combining advanced technologies in one 5-layer dressing, ALLEVYN COMPLETE CARE Dressing helps to protect fragile⁵ skin and minimise disruption^{3,4} – optimising patient comfort⁶ and restoring confidence in care.

UNIQUE NEEDS, COMPLETE CARE



What does ALLEVYN[◇] COMPLETE CARE Foam Dressing mean for you and your team?

It means caring – with confidence

- ✓ **Confidence** that your dressing will stay in place^{6,7}
- ✓ **Confidence** in achieving longer wear times up to 7 days*⁸
- ✓ **Confidence** in protecting fragile skin⁵
- ✓ **Confidence** in undisturbed healing progression^{3,4}

With fewer unplanned dressing changes and fewer interruptions, you can have time back to focus on what really matters.

And as for your patients, they get the freedom to live the life that matters to them.

Patients deserve a life with improved comfort, less anxiety and fear of leakage.[†]

Less focus on dressing changes,^{6,7} and more focus on helping patients – this is wound management as it should be.

Reclaim your time and focus on what truly counts

*As demonstrated in a volunteer trial. 120 volunteers wore the dressing on intact skin for up to 7 days on their thighs or shins
†Not clinically proven



Patients like Zoe can enjoy getting back

to the life they love

Every ALLEVYN[®] COMPLETE CARE Foam Dressing contains enhanced technologies to **minimise strike-through**,⁹ enable extended wear times of up to 7 days^{6,7} and support undisturbed healing.^{3,4}



Visually masks unsightly exudate to keep the dressing **looking fresh for longer** and shows when it needs changing.^{8,10}



No need to change



Consider changing



Time to change



ExuMASK[®]
Technology

Unbonded design enables layers to frictionally slide within the dressing to **absorb shear forces and which may help reduce soft tissue strain**^{*11} and protects against friction, pressure, and microclimate.^{†12-14}



ShearDEFENSE[®]
Technology



ExuLOCK[®]
Technology

Provides superior absorbency^{||15,16} that **locks in exudate and bacteria**,^{17,18} maintains an optimal wound healing environment,^{‡12,19-21} and works effectively without primarily relying on evaporation.^{21,22}



HighFLEX[®]
Technology

Multi-way stretch **helps to fit around awkward body contours**.^{19,23}



ComfortSTAY[®]
Technology

Our next-gen adhesive provides secure hold with gentle adhesion^{3,5,24,25} and **allows easy repositioning for skin checks**.^{24§*}

*Tested in-vitro †As demonstrated in silico and in vitro ‡Over the wear time of the dressing §p<0.001; As demonstrated in vitro ||compared to leading competitors

Ordering information:

Shape	S+N Code	Size	Description	Pad size	Carton
	66807695	7.5 x 7.5cm	Square XS	4.5 x 4.5cm	10
	66807696	10 x 10cm	Square S	6.5 x 6.5cm	10
	66807697	12.5 x 12.5cm	Square M	8.5 x 8.5cm	10
	66807698	15 x 15cm	Square L	11 x 11cm	10
	66807699	17.5 x 17.5cm	Square XL	13.5 x 13.5cm	10
	66807705	25 x 25.2cm	Heel	20 x 20.2cm	10
	66807703	17.2 x 17.4cm	Sacrum S	8.4 x 12.4cm	10
	66807704	21.5 x 23cm	Sacrum L	12.3 x 16.9cm	10
	66807700	10 x 20cm	Rectangle S	5 x 15cm	10
	66807701	10 x 30cm	Rectangle L	5 x 25cm	10
	66807702	17.1 x 17.8cm	Multisite	11.5 x 12.3cm	10

To order or
find out more,
scan this code:



References: **1.** Atkin L. Clinical inertia in chronic wound care. Wounds International. 2025;16(1):48-51. **2.** Li Z, Lin F, Thalib L, Chaboyer W. Global prevalence and incidence of pressure injuries in hospitalised adult patients: a systematic review and meta-analysis. Int J Nurs Stud. 2020;105:103546. **3.** Vowden K, Moiem N, Drysdale K, Mistry C. An open, prospective randomised, multi-centre clinical evaluation of a hydrocellular silicone foam dressing* in the management of exuding chronic and acute wounds. Paper presented at: EWMA; 2011; Brussels. **4.** Hurd T, Gregory L, Jones A, Brown S. A multi-centre in-market evaluation of ALLEVYN® Gentle Border. Wounds UK. 2009;5(3):32 - 44. **5.** Smith+Nephew. Internal Report: EA/AWM/GENERAL/001/v2. **6.** Smith+Nephew 2025. Topline results from ALLEVYN® COMPLETE CARE Foam Dressing volunteer trial (HVS2313) on thighs and shins. Internal report: CSD.AWM.25.032. **7.** Atkinson L, Allen D, Costa B. Poster presented at: European Wound Management Association (EWMA); May 1–3, 2024; London, UK. (NB: EM meta-analysis) **8.** Costa B, Allen D. Results from an international survey: foam dressings with WTCI may enhance clinician confidence to extend wear times compared with other foam dressings. Poster presented at: the 33rd conference of the EWMA; May 3-5, 2023; Milan, Italy. **9.** Smith+Nephew 2025. Wound Models testing over 4 Days – Horizontal, Moderate Flow, Viscous SWF and Wound Models Testing over 7 Days, Horizontal, Moderate Flow of ALLEVYN COMPLETE CARE Dressing. Internal report: CSD.AWM.25.009 V2. **10.** Smith+Nephew 2025. Change Indicator Usability of ALLEVYN COMPLETE CARE. Internal report: CSD.AWM.25.023 V2. **11.** Orlova D, Orlov A, Gefen A. The Protective Efficacy of a New Soft Silicone Multi-Layer Dressing in Reducing the Heel Pressure Ulcer Risk. Int Wound J. 2025;22(10):e70764. **12.** Smith+Nephew 2025. Internal report: CSD.AWM.25.011. **13.** Smith+Nephew 2024. Internal report: CSD.AWM.24.057 V2. **14.** Smith+Nephew 2025. Internal report: CSD.AWM.25.016 **15.** Smith+Nephew 2025. Internal report: CSD.AWM.25.079 **16.** Smith+Nephew 2026. Internal report: CSD.AWM.26.001. **17.** Smith+Nephew 2025. Free swell Absorbency and Retention after Compression Testing of ALLEVYN COMPLETE CARE Dressing. Internal report: CSD.AWM.25.008. **18.** Smith+Nephew 2025. Bacterial absorbance, sequestration & retention testing of ALLEVYN COMPLETE CARE. Internal report: CSD.AWM.25.019 V2. **19.** Rossington A, Drysdale K, Winter R. Clinical performance and positive impact on patient wellbeing of ALLEVYN Life. Wounds UK. 2013;9(4):91 - 95. **20.** Smith+Nephew 2025. Internal report: CSD.AWM.25.009 V2. **21.** Smith+Nephew 2025. Internal report: CSD.AWM.25.010. **22.** Smith+Nephew 2025. Internal report: CSD.AWM.25.008. **23.** Stephen-Haynes J, Bielby A, Searle R. The clinical performance of a silicone foam in an NHS community trust. JCN. 2013;27(5):50-59. **24.** Smith+Nephew 2025. Internal report: CSD.AWM.25.015 V2. **25.** Smith+Nephew 2025. Internal report: CSD.AWM.25.059.

Advanced Wound Management
Smith & Nephew, Medical Ltd.
101 Hesse Road
Hull HU3 2BN, UK
T +44 (0) 1482 225181
F +44 (0) 1482 328326

www.smith-nephew.com

© Trademark of Smith+Nephew
All Trademarks acknowledged
© March 2026 Smith+Nephew
AWM-AWC-50459 | GMC2037

For detailed product information, including indications for use, contraindications, warnings and precautions, please consult the product's Instructions for Use (IFU) prior to use. Products may not be available in all markets because product availability is subject to the regulatory and / or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area.