



Tray Layout




TRIGEN Base Instrument Set

Set No. 7167-4012

Cat. Item	Description	Qty	Reference
7112-9401	Small Outer Case	1	N/A
7112-9402	Lid for Outer Case	1	N/A
7167-4021	TRIGEN Base Tray	1	N/A
7163-1066	Medium Hexdriver	1	A
7163-1068	Short Hexdriver	1	B
7163-1116	12.5mm Entry Reamer	1	C
7163-1140	Guide Bolt Wrench	1	D
7163-1152	9.0mm Drill Sleeve	2	E
7163-1161	Multipurpose Driver	1	F
7163-1186	Mini Connector	1	G
7163-1189	Screw Depth Gauge	1	H
7167-4000	Cannulated Awl	1	I
7167-4060	Entry Portal Tube	1	J
7167-4074	3.2mm T-handle Trocar	1	K
7167-4075	Honeycomb	1	L
7167-4076 or 7167-4576	T-handle	1	M
7167-4077	Reducer	1	N
7167-4077	Reducer Connector	1	O

Cat. Item	Description	Qty	Reference
7167-4078	Obturator	1	P
7167-4079	Ruler	1	Q
7167-4080	Gripper	1	R
7167-4081	Impactor	1	S
7167-4082	Slotted Hammer	1	T
7167-4083	4.0mm Drill Sleeve	2	U
7167-4084	Screwdriver Release Handle	1	V
7167-4085	Screw Length Sleeve	1	W
7167-4092	Entry Portal Handle	1	X
7167-1212	TRIGEN Reamer Set, Optional	1	Y
7111-8200	SCULPTOR Flexible Reamer, Optional	1	Z
7163-1130	Flexible Reamer Extender, Optional	1	AA
7164-1123	TRIGEN 4.0 Diaphyseal Drill	1	BB
7163-1121	4.0mm Long AO Pilot Drill, 333mm, Disposable	2	CC
7163-1123	4.00mm Short AO Pilot Drill, 161mm, Disposable	1	DD
7167-4130	3.2mm x 343mm Brad Point Tip Guide Pin, Disposable	3	EE
7163-1070	Long Hexdriver, Optional	1	FF
7175-1153	AO Mini Connector, Optional	1	GG
7163-1187	Trinkle to Mini Connector, Optional	1	HH



TRIGEN[®] Percutaneous TAN[®] FAN Instrument Set
 Set No. 7163-2351

Cat. Item	Description	Qty	Reference
7163-1021	Percutaneous Drill Guide	1	A
7163-1022	Radiolucent Drop	1	B
7163-1039	14mm Channel Reamer	1	C
7163-1024	Percutaneous Guide Bolt	2	D
7163-1036	Sterile 6.4mm Tap	1	E
7163-1026	4.0mm Trocar Drill Sleeve	2	F
7163-1191	4.0mm Trocar	1	G
7163-1025	AP Alignment Tower	1	H
7163-1015	AP Alignment Arm	1	I
7163-1035	Sterile 6.4mm Drill	1	J
7163-1320	TRIGEN Disposable Nail Extractor	1	K

Smith & Nephew, Inc.
 1450 Brooks Road
 Memphis, TN 38116
 USA

www.smith-nephew.com

Telephone: 1-901-396-2121
 Information: 1-800-821-5700
 Orders and Inquiries: 1-800-238-7538

[®]Trademark of Smith & Nephew. Certain marks Reg. US Pat. & TM Off.
 All Trademarks acknowledged.

©2015 Smith & Nephew.
 00283 V3 04/15