



New!

Get on the right track
with first-time access


Combine with our new ARTHROGARDE® Hip Access Cannulas for optimal protection.

 **smith&nephew**
CROSSTRAC◊
Hip Guide System



Achieve accurate trajectories and first-time access to the hip joint using our new CROSSTRAC[®] Hip Guide System*

Accessing the hip joint has never been easier. With its unique triangulated design and needle point accuracy, the CROSSTRAC Hip Guide System allows for consistent trajectory alignment and portal placement where you want it. Use with our standard length 4mm arthroscopes and our new atraumatic-in-design ARTHROGARDE[®] Hip Access Cannulas for improved protection and safety.




Scope – based design

Specialized connection hub between scope/cannula assembly and the guide assures repeatable triangulation under direct visualization



Accurate Trajectories

Aligns to any portal position within its range of motion and cannulated bullet assures accurate trajectories are achieved



First-time Access

Decrease the number of attempts to establish portal placement and gain access where you want it

Ordering Information

72201832 CROSSTRAC Hip Guide System

System Includes:

72201831 CROSSTRAC Hip Guide with Bullet
72201741 4.5 mm ARTHROGARDE Hip Access Cannula GREEN (CROSSTRAC Hub)
72201742 5.0 mm ARTHROGARDE Hip Access Cannula BLUE (CROSSTRAC Hub)
72200096 4.5 mm Cannulated Obturator with Delrin Handle
72200097 5.0 mm Cannulated Obturator with Delrin Handle

Sold Separately:

72201811 XL (extra-long) Disposable Arthroscopy Needle



Endoscopy
Smith & Nephew, Inc.
Andover, MA 01810
USA

www.smith-nephew.com
T +1 978 749 1000
F +1 978 749 1108
US Customer Service: +1 800 343 5717

*Trademarks of Smith & Nephew.
Certain marks Reg US Pat & TM off.
©2008 Smith & Nephew, Inc. All rights reserved.
Printed in USA. 04/08 10600396 Rev A

*Patent Pending