1. Identification

1.1. Product identifier

Product Identity

Acticoat* (with SILCRYST Nanocrystals) Absorbent Dressing * Trademark of Smith & Nephew - ™ SILCRYST is a trademark of Smith & Nephew , Inc.

Alternate Names Acticoat Absorbent Dressing with silver coating

1.2. Relevant identified uses of the substance or mixture and uses advised against

Medical device used for infection control. Intended use

Application Method See Technical Data Sheet.

1.3. Details of the supplier of the safety data sheet

Company Name Smith & Nephew

970 Lake Carillon Drive, Suite 110

St. Petersburg, FL 33716

Emergency

Customer Service: Smith & Nephew 1-800-876-1261

2. Hazard(s) identification

2.1. Classification of the substance or mixture

Eye Dam. 1;H318 Causes serious eye damage.

Aquatic Chronic 1;H410 Very toxic to aquatic life with long lasting effects.

2.2. Label elements

Using the Toxicity Data listed in section 11 and 12 the product is labeled as follows.



Danger

H318 Causes serious eye damage.

H410 Very toxic to aquatic life with long lasting effects.

[Prevention]:

P273 Avoid release to the environment.

P280 Wear protective gloves / eye protection / face protection.

[Response]:

P305+351+338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do - continue rinsing.

P310 Immediately call a POISON CENTER or doctor / physician.

P391 Collect spillage.

[Storage]:

No GHS storage statements

[Disposal]:

P501 Dispose of contents / container in accordance with local / national regulations.

3. Composition/information on ingredients

This product contains the following substances that present a hazard within the meaning of the relevant State and Federal Hazardous Substances regulations.

Ingredient/Chemical Designations	Weight %	GHS Classification	Notes
Alginic acid, calcium salt CAS Number: 0009005-35-0	75 - 100	Not Classified	[1]
Silver CAS Number: 0007440-22-4	1.0 - 10	Not Classified	[1][2]
Silver oxide CAS Number: 0020667-12-3	1.0 - 10	Ox. Sol. 1;H271 Eye Dam. 1;H318 Aquatic Acute 1;H400 Aquatic Chronic 1;H410	[1]

In accordance with paragraph (i) of §1910.1200, the specific chemical identity and/or exact percentage (concentration) of composition has been withheld as a trade secret.

4. First aid measures

4.1. Description of first aid measures

General In all cases of doubt, or when symptoms persist, seek medical attention.

Never give anything by mouth to an unconscious person.

Inhalation This dressing does not present an inhalation hazard.

Eyes Irrigate copiously with clean water for at least 15 minutes, holding the eyelids apart and

seek medical attention.

Skin Obtain medical attention in the event skin irritation develops or persists. **Ingestion** Ingestion of the dressing will act as a foreign body. Obtain medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Overview <u>Acute Effects</u>

Temporary discoloration of affected areas or rashes and irritation in sensitized individuals.

Eye: Corneal irritation may occur as a result of mechanical abrasion. In the event of eye

irritation, flush eyes with flowing water.

Skin: Skin irritation may result as an allergic reaction to the silver content of the dressing.

Obtain medical attention in the event irritation develops or persists.

Ingestion: Ingestion of the dressing may cause a choking hazard or internal blockage.

Obtain medical attention.

Chronic Effects

^[1] Substance classified with a health or environmental hazard.

^[2] Substance with a workplace exposure limit.

^[3] PBT-substance or vPvB-substance.

^{*}The full texts of the phrases are shown in Section 16.

Excessive or chronic exposure to silver by skin contact, ingestion or inhalation may cause Argyrosis - a local or generalized discoloration of the eyes, skin and mucous membranes, where grey or blue patches of pigmentation are formed.

See section 2 for further details.

Eyes Causes serious eye damage.

5. Fire-fighting measures

5.1. Extinguishing media

Dressing is not reactive to water or other means of extinguishment. If involved in a fire, treat as a flammable solid.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition: No hazardous decomposition data available.

5.3. Advice for fire-fighters

A Self-Contained Breathing Apparatus should be worn.

Do not release runoff from fire control methods to sewers or waterways.

ERG Guide No. ----

6. Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Put on appropriate personal protective equipment (see section 8).

6.2. Environmental precautions

Do not allow spills to enter drains or waterways.

Use good personal hygiene practices. Wash hands before eating, drinking, smoking or using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

6.3. Methods and material for containment and cleaning up

Treat used dressing as medical waste. No special spill procedures required.

7. Handling and storage

7.1. Precautions for safe handling

No special handling procedures required. If used, handle as per medical waste.

See section 2 for further details. - [Prevention]:

7.2. Conditions for safe storage, including any incompatibilities

Handle containers carefully to prevent damage and spillage.

Incompatible materials: The silver coating is incompatible with these products: Acetylene, Acids, Amines, Ammonia, Hydrazine, Metal Sulfides, Peroxides, Strong oxidizers.

Store in a cool, dry environment. Keep away from moisture. Avoid extreme heat and strong oxidizing agents.

See section 2 for further details. - [Storage]:

7.3. Specific end use(s)

No data available.

8. Exposure controls and personal protection

8.1. Control parameters

Exposure

CAS No.	Ingredient	Source	Value
0007440-22-4	Silver	OSHA	TWA 0.01 mg/m³
		ACGIH	TWA: 0.01 mg/m3STEL: 0.03 mg/m ³
		NIOSH	TWA 0.01 mg/m ³
		Supplier	No Established Limit
0009005-35-0 Alginic acid, calcium salt	OSHA	No Established Limit	
		ACGIH	No Established Limit
	NIOSH	No Established Limit	
	Supplier	No Established Limit	
0020667-12-3	Silver oxide	OSHA	No Established Limit
		ACGIH	No Established Limit
		NIOSH	No Established Limit
		Supplier	No Established Limit

Carcinogen Data

CAS No.	Ingredient	Source	Value
0007440-22-4	0007440-22-4 Silver		Select Carcinogen: No
		NTP	Known: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;
0009005-35-0 Alginic acid, calcium salt		OSHA	Select Carcinogen: No
		NTP	Known: No; Suspected: No
	IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;	
0020667-12-3	Silver oxide		Select Carcinogen: No
		NTP	Known: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;

8.2. Exposure controls

RespiratoryNo special protection is required.EyesNo special protection is required.SkinNo special protection is required.

Engineering Controls None required.

Other Work Practices Use good personal hygiene practices. Wash hands before eating, drinking, smoking or

using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

See section 2 for further details. - [Prevention]:

9. Physical and chemical properties

Appearance Grey/Brown Non-woven Fabric

Odor Odorless

Odor thresholdNot MeasuredpHNot availableMelting point / freezing pointNot applicableInitial boiling point and boiling rangeNot applicable

Flash Point

Evaporation rate (Ether = 1)

Flammability (solid, gas)

Not applicable

Not Applicable

Upper/lower flammability or explosive limits Lower Explosive Limit: Not available

Upper Explosive Limit: Not available

Vapor pressure (Pa)Not applicableVapor DensityNot applicable

Specific Gravity > 1

Solubility in Water Silver coating is sparingly soluble in water.

Auto-ignition temperatureNot availableDecomposition temperatureNot availableViscosity (cSt)Not applicableVOC ContentNot availableTasteMetallic

9.2. Other information

No other relevant information.

10. Stability and reactivity

10.1. Reactivity

Hazardous Polymerization will not occur.

10.2. Chemical stability

Stable under normal circumstances.

10.3. Possibility of hazardous reactions

Fiber will gel or dissolve in an aqueous solution containing sodium or certain ions.

10.4. Conditions to avoid

Avoid extreme heat and strong oxidizing agents.

10.5. Incompatible materials

The silver coating is incompatible with these products: Acetylene, Acids, Amines, Ammonia, Hydrazine, Metal Sulfides, Peroxides, Strong oxidizers.

10.6. Hazardous decomposition products

No hazardous decomposition data available.

11. Toxicological information

Acute toxicity

Ingredient	Oral LD50, mg/kg	Skin LD50, mg/kg	Inhalation Vapor LC50, mg/L/4hr	Inhalation Dust/Mist LC50, mg/L/4hr	Inhalation Gas LC50, ppm
Alginic acid, calcium salt - (9005-35-0)	No data available	No data available	No data available	No data available	No data available
Silver - (7440-22-4)	No data available	No data available	No data available	No data available	No data available
Silver oxide - (20667-12-3)	No data available	No data available	No data available	No data available	No data available

Note: When no route specific LD50 data is available for an acute toxin, the converted acute toxicity point estimate was used in the calculation of the product's ATE (Acute Toxicity Estimate).

Classification	Category	Hazard Description
Acute toxicity (oral)		Not Applicable
Acute toxicity (dermal)		Not Applicable
Acute toxicity (inhalation)		Not Applicable
Skin corrosion/irritation		Not Applicable
Serious eye damage/irritation	1	Causes serious eye damage.
Respiratory sensitization		Not Applicable
Skin sensitization		Not Applicable
Germ cell mutagenicity		Not Applicable
Carcinogenicity		Not Applicable
Reproductive toxicity		Not Applicable
STOT-single exposure		Not Applicable
STOT-repeated exposure		Not Applicable
Aspiration hazard		Not Applicable

12. Ecological information

12.1. Toxicity

Very toxic to aquatic life with long lasting effects.

Toxic to aquatic life

Aquatic Ecotoxicity

Ingredient	96 hr LC50 fish, mg/l	48 hr EC50 crustacea, mg/l	Er C50 algae, mg/l
Alginic acid, calcium salt - (9005-35-0)	Not Available	Not Available	Not Available
Silver - (7440-22-4)	Not Available	Not Available	Not Available
Silver oxide - (20667-12-3)	Not Available	Not Available	Not Available

12.2. Persistence and degradability

There is no data available on the preparation itself.

12.3. Bioaccumulative potential

Not Measured

12.4. Mobility in soil

No data available.

12.5. Results of PBT and vPvB assessment

This product contains no PBT/vPvB chemicals.

12.6. Other adverse effects

No data available.

13. Disposal considerations

13.1. Waste treatment methods

Observe all federal, state and local regulations when disposing of this substance.

14. Transport information

DOT (Domestic Surface IMO / IMDG (Ocean ICAO/IATA Transportation) Transportation) 14.1. UN number Not Applicable Not Regulated Not Regulated 14.2. UN proper shipping Not Regulated Not Regulated Not Regulated name 14.3. Transport hazard **DOT Hazard Class:** Not **IMDG:** Not Applicable Air Class: Not Applicable class(es) Applicable Sub Class: Not Applicable 14.4. Packing group Not Applicable Not Applicable Not Applicable

14.5. Environmental hazards

IMDG Marine Pollutant: Yes (Silver oxide)

14.6. Special precautions for user

No further information

15. Regulatory information

Regulatory Overview The regulatory data in Section 15 is not intended to be all-inclusive, only selected

regulations are represented.

Toxic Substance All components of this material are either listed or exempt from listing on the TSCA

Control Act (TSCA) Inventory.

WHMIS Classification D2B E

US EPA Tier II Hazards Fire: No

Sudden Release of Pressure: No

Reactive: No

Immediate (Acute): Yes

Delayed (Chronic): No

EPCRA 311/312 Chemicals and RQs (lbs):

Silver (1,000.00)

EPCRA 302 Extremely Hazardous:

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

EPCRA 313 Toxic Chemicals:

Silver

Silver oxide

Proposition 65 - Carcinogens (>0.0%):

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

Proposition 65 - Developmental Toxins (>0.0%):

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

Proposition 65 - Female Repro Toxins (>0.0%):

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

Proposition 65 - Male Repro Toxins (>0.0%):

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

New Jersey RTK Substances (>1%):

Silver

Pennsylvania RTK Substances (>1%):

Silver

16. Other information

The information and recommendations contained herein are based upon data believed to be correct. However, no guarantee or warranty of any kind, expressed or implied, is made with respect to the information contained herein. We accept no responsibility and disclaim all liability for any harmful effects which may be caused by exposure to our products. Customers/users of this product must comply with all applicable health and safety laws, regulations, and orders.

The full text of the phrases appearing in section 3 is:

H271 May cause fire or explosion; strong oxidizer.

H318 Causes serious eye damage.

H400 Very toxic to aquatic life.

H410 Very toxic to aquatic life with long lasting effects.

Not Classified Not Classified

This is the first version in the GHS SDS format. Listings of changes from previous versions in other formats are not applicable.

End of Document