

Resection and Fluid Management Technologies

DYONICS[®] Blades and Burrs | COBLATION[®] Wands

Smith+Nephew

MECHANICAL RESECTION

PLATINUM BLADES		Size	Ref #	Hub Color	Max Fwd/Rev Speed*
	PLATINUM FLYER [®] Blade	4.0mm	72205292		8,000**
	BONECUTTER [®] PLATINUM Blade	3.5mm 4.5mm 5.5mm	72205108 72202531 72202530	  	10,000** 10,000** 10,000**
	INCISOR PLUS PLATINUM Blade	3.5mm 4.5mm 5.5mm	72205106 72203013 72203519	  	5,000 5,000 5,000
	Curved INCISOR [®] Plus PLATINUM Blade	4.5mm	72205109		5,000
	SYNOVATOR [®] PLATINUM Blade	3.5mm 4.5mm 5.5mm	72205107 72203523 72203524	  	5,000 5,000 5,000
	Curved SYNOVATOR PLATINUM Blade	4.5mm	72205110		5,000
	PLATINUM FR Blade	4.5mm 5.5mm	72203521 72203522	 	5,000 5,000

ESSENTIAL BLADES		Size	Ref #	Hub Color	Max Fwd/Rev Speed*
	BONECUTTER Blade (Full Radius)	4.5mm 5.5mm	7206011 7206010	 	8,000** 8,000**
	INCISOR PLUS ELITE Blade	3.5mm 4.5mm 4.5mm CONCV 5.5mm	72200095 7210976 72200036 72200081	   	5,000 5,000 5,000 5,000
	FULL RADIUS Blade	3.5mm 3.5mm CONCV 4.5mm 4.5mm ORBIT 5.5mm	7205305 7210980 7205306 7205323 7205307	    	5,000 3,000 5,000 5,000 5,000
	FULL RADIUS ELITE [®] Blade	3.5 mm 4.5mm 4.5mm CONCV 5.5mm	7210751 7210499 7210988 7210909	   	5,000 5,000 5,000 5,000
	RAZORCUT Blade	3.5mm 4.5mm 5.5mm	7205317 7205318 7205319	  	5,000 5,000 5,000
	SYNOVATOR Blade	4.5mm 4.5mm CONCV 4.5mm ORBIT 5.5mm	7205310 7205334 7205321 7205311	   	5,000 5,000 5,000 5,000
	CUTTER Blade	3.5mm 4.5mm	7205308 7205309	 	5,000 5,000
	TURBOWHISKER Blade	4.5mm	7205316		3,000
	ELECTROBLADE [®] ELITE Blade (INCISOR)	4.5mm	7209983		5,000
	ELECTROBLADE ELITE Blade (INCISOR)	4.5mm	7209700		5,000
	ELECTROBLADE BONECUTTER Blade (Full Radius)	4.5mm 5.5mm	72202213 7205962	 	8,000** 8,000**
	ELECTROBLADE BONECUTTER Blade (Full Radius)	4.5mm 4.5mm	7209855 7205961	 	5,000 5,000

STANDARD BURRS		Size	Ref #	Hub Color	Max Fwd/Rev Speed*
	HELICUT [®] Burr	4.0mm	7205727		5,000
	ELITE ABRADER Burr 10k RPM	4.0mm 5.5mm	72200722 72200723	 	10,000 10,000
	ELITE NOTCHBLASTER [®] Burr 10k RPM	4.0mm 5.5mm	72200726 72200727	 	10,000 10,000
	ELITE ACROMIONIZER Burr 10k RPM	4.0mm 5.5mm	72200724 72200725	 	10,000 10,000
	ELITE ACROMIOBLASTER [®] Burr 10k RPM	4.0mm 5.5mm	72200730 72200731	 	10,000 10,000
	ELITE STONECUTTER Burr 10k RPM	4.0mm 5.5mm	72200728 72200729	 	10,000 10,000

HIP BLADES		Size	Ref #	Hub Color	Max Fwd/Rev Speed*
	BONECUTTER Blade (Full Radius)	4.5mm	7210979		8,000**
	INCISOR PLUS ELITE Blade	4.5mm 4.5mm CONCV	72200414 72200494	 	3,000 3,000
	FULL RADIUS Blade	4.5mm 4.5mm CONCV	7205341 7205335	 	3,000 3,000

HIP BURRS		Size	Ref #	Hub Color	Max Fwd/Rev Speed*
	ABRADER Burr	4.0mm 5.5mm	72200080 72200082	 	8,000 8,000
	ABRADER Burr High Visibility Sheath	4.0mm 5.5mm	72203128 72203127	 	8,000 8,000
	BARREL Burr High Visibility Sheath	4.0mm 5.5mm	72203132 72203131	 	8,000 8,000
	FLAT TOP Burr High Visibility Sheath	4.0mm 5.5mm	72203130 72203129	 	8,000 8,000


EXTREMITIES BLADES		Size	Ref #	Hub Color	Max Fwd/Rev Speed*
	CUTTER Blade	2.9mm	72201506		5,000**
	FULL RADIUS Blade	2.0mm 2.9mm 3.5mm	72201507 72201509 72201510	  	5,000** 5,000** 5,000**
	INCISOR PLUS ELITE Blade	2.9mm 3.5mm	72201513 72201514	 	5,000** 5,000**

EXTREMITIES BURRS		Size	Ref #	Hub Color	Max Fwd/Rev Speed*
	ABRADER Burr	2.9mm 3.5mm	72201518 72201519	 	6,000 6,000
	BARREL Burr	2.9mm	72201520		6,000

*Indicates the max forward and reverse speed for increased resection rate. Each blade operates at a speed of 8 unless otherwise indicated while in the proprietary oscillate mode. In this mode, the blade revolves two times in each direction during an oscillate period. The window opens frequently, drawing and rapidly resecting soft tissue.

**Operates at a speed of 9 in proprietary oscillate mode for increased resection rate.

COBLATION

SHOULDER		Ref #
	FLOW 90 [°] Wand	72290038 <input type="checkbox"/>
	AMBIENT [®] SUPER TURBOVAC [®] 90 IFS Wand	ASHA4250-01 <input type="checkbox"/> <input type="checkbox"/>
	SUPER TURBOVAC 90 IFS Wand	ASH4250-01 <input type="checkbox"/> <input type="checkbox"/>
	SUPER TURBOVAC 90 ICW Wand	ASC4250-01 <input type="checkbox"/> <input type="checkbox"/>
	AMBIENT MEGAVAC 90 ICW Wand	ASCA5001-01 <input type="checkbox"/> <input type="checkbox"/>
	STARVAC [®] 90 ICW Wand	ASC4251-01 <input type="checkbox"/> <input type="checkbox"/>
	COVATOR [®] 20 Wand	AC4340-01 <input type="checkbox"/> <input type="checkbox"/>
	TOPAZ [®] XL IFS Wand	ACH4045-01 <input type="checkbox"/> <input type="checkbox"/>

KNEE		Ref #
	FLOW 50 [°] Wand	72290037 <input type="checkbox"/>
	AMBIENT [®] SUPER MULTIVAC [®] 50 IFS Wand	ASHA4830-01 <input type="checkbox"/> <input type="checkbox"/>
	SUPER MULTIVAC 50 IFS Wand	ASH4830-01 <input type="checkbox"/> <input type="checkbox"/>
	SUPER MULTIVAC 50 ICW Wand	ASC4830-01 <input type="checkbox"/> <input type="checkbox"/>
	AMBIENT COVAC [®] 50 IFS Wand	ASHA2530-01 <input type="checkbox"/> <input type="checkbox"/>
	AMBIENT COVAC 70 IFS Wand	ASHA3730-01 <input type="checkbox"/> <input type="checkbox"/>
	PARAGON T2 [®] Wand	AC5331-01 <input type="checkbox"/>
	SABER [®] 30 Wand	AC4330-01 <input type="checkbox"/>

HIP		Ref #
	AMBIENT HIPVAC [®] 50 IFS Wand	72290004 <input type="checkbox"/> <input type="checkbox"/>
	MULTIVAC 50 XL ICW Wand	ASC4730-01 <input type="checkbox"/>
	TURBOVAC 90 XL ICW Wand	ASC1336-01 <input type="checkbox"/>

EXTREMITIES		Ref #
	Short BEVEL 35 DEGREE ICW Wand	AC2823-01 <input type="checkbox"/>
	MICROBLATOR [®] 30 ICW Wand	AC4050-01 <input type="checkbox"/>
	TOPAZ [®] EZ IFS Wand	ACH4041-01 <input type="checkbox"/> <input type="checkbox"/>

IFS = finger switches, ICW = no finger switches

EQUIPMENT

	WEREWOLF ⁺ COBLATION [®] System	72290146
	WEREWOLF Wired Footswitch	72290007
	WEREWOLF Wireless Footswitch	72290008
	QUANTUM [®] 2 Controller (H4500-00)	28168
	QUANTUM 2 ATLAS Wired Foot Switch	10863
	DYONICS [®] POWER II Control Unit	72200873
	DYONICS POWERMAX ELITE Handpiece	72200616
	DYONICS PLATINUM Handpiece	72205321
	DYONICS POWERMINI Handpiece	72201500
	DYONICS POWER II Footswitch	72201092
	DOUBLEFLO [®] Pump and Manual [†]	72205776 (Global Except EU) 72205775 (EU Only)
	Day Tube Set [†] , 15 per box	72205772
	Patient Tube Set [†] , 20 per box	72205773
	Inflow Tube Set [†] , 15 per box	72205774
	Outflow Tube Set [†] , 15 per box	72205765

Compatible with:

- WEREWOLF and WEREWOLF+ COBLATION Systems
- QUANTUM 2 System

Learn more at smith-nephew.com

Distributed by:
Sports Medicine
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810

www.smith-nephew.com
T +978 749 1000
US Customer Service:
+1 800 343 5717

Manufactured by:
ArthroCare Corporation
7000 West William Cannon Dr.
Austin, TX 78735, USA

Manufactured by:

85 rue du Chêne vert
31670 Labège – France

©Trademark of Smith+Nephew.
©2024 Smith & Nephew. All rights reserved. All trademarks acknowledged. Printed in USA. 09809 V4 05/24