

Maximum flexibility + minimal effort

A new standard for controlled ankle joint distraction

Smith+Nephew

SPIDER2[◇]
Limb Positioner



Controlled distraction

With the comfort of personalized ankle positioning

The SPIDER2[◇] Ankle set-up allows the user to move the ankle in their most preferred position to ensure appropriate access during their procedure. Giving the user the ability to re-position and adjust traction with ease intraoperatively, allows the SPIDER2 Ankle set-up to bring a new dimension to ankle procedures.



KEY FEATURES

Maintaining a steady position from start to finish, the SPIDER2 Limb Positioner is designed to allow surgeons full control of limb positioning, range of motion and traction from within the sterile field.¹

SPIDER2 Limb Positioner

- Provides quick and simple set-up²⁻⁴
- Designed to control knee flexion angle to assist with appropriate ankle positioning²⁻⁴



Traction Accessory

- Indicator for increasing and decreasing fine traction
- Indicator for measurement of traction
- Simple process for achieving controlled traction



Ankle Arthroscopy Essentials

- Locking mechanism for disposable strap provides tactile click to know when it is locked into position
- Sterile drape allows for SPIDER2 Limb Positioner to be manipulated by the user intraoperatively with ease¹



FERKEL™ Thigh Holder

- Provides the necessary counter traction needed for foot and ankle arthroscopy procedures while also stabilizing the thigh
- Designed to clamp to the bedrail of a surgical table, the FERKLE Thigh Holder comes with disposable pads for easy clean-up



Ordering information

SPIDER2° Ankle Distraction System	
Reference	Description
72203299	SPIDER2 Limb Positioner
7210577	Traction Accessory
7210570	Piggy Back Connector, box of 2
72202932	Ankle Arthroscopy Accessory, 2 units
72201812	FERKEL™ Thigh Holder
Consumables	
Reference	Description
72202931	Ankle Stabilization Kit, box of 6
72201813	FERKEL™ Thigh Holder, box of 1

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Learn more at smith-nephew.com

Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

www.smith-nephew.com
T +978 749 1000
US Customer Service:
+1 800 343 5717

°Trademark of Smith+Nephew.
All trademarks acknowledged.
©2024 Smith+Nephew. All rights reserved. 43153 V1 09/24

References

1. Tenet Medical Engineering 2013. Design Verification Report - Distal Thigh Holder Rev. B. 2. Tenet Medical Engineering 2008. Design Verification Results - Spider Leg Accessory Rev A. 3. Tenet Medical Engineering 2009. Design Validation Protocol - Ankle Arthroscopy Accessory and Ankle Stabilization Kit 4. Tenet Medical Engineering 2011 Design Verification Results Ankle Arthroscopy Accessory Revision B