
Product ordering codes

Smith+Nephew

Dressings	Size	Fluid handling capacity[1]	PICO◊ 7 device + 1 dressing Code	PICO◊ 7 device + 2 dressings Code	PICO 14 device + 2 dressings Code	Multipack with 5 dressings Code	PICO 7Y device + 2 dressings Code
	Multisite small 15cm x 20cm	Up to 135ml*	66802010	66802000	66802040	66802020	–
	Multisite large 20cm x 25cm	Up to 272ml*	66802011	66802001	66802041	66802021	66802031
	9cm x 15cm	Up to 44ml*	–	66807749	66807754	66807757	–
	10cm x 20cm	Up to 82ml*	66802012	66802002	66802042	66802022	–
	10cm x 25cm	Up to 110ml*	66807752	66807750	66807755	66807758	–
	10cm x 30cm	Up to 137ml*	66802013	66802003	66802043	66802023	–
	10cm x 40cm	Up to 192ml*	66802014	66802004	66802044	66802024	–
	10cm x 55cm	Up to 275ml*	66807753	66807751	66807756	66807759	–
	15cm x 15cm	Up to 110ml*	66802015	66802005	66802045	66802025	–
	15cm x 20cm	Up to 165ml*	66802016	66802006	66802046	66802026	–
	15cm x 30cm	Up to 275ml*	66802017	66802007	66802047	66802027	–
	20cm x 20cm	Up to 247ml*	66802018	66802008	66802048	66802028	–
	25cm x 25cm	Up to 440ml*	66802019	66802009	66802049	66802029	–

Consumables	Size	Code
 Foam dressing filler	10cm x 12.5cm	66801021
 5 Antimicrobial Gauze Rolls + 1 SECURA◊ NSBF Wipe	11.4cm x 3.7m	66802127

*Derived from a bench testing of a PICO sNPWT Dressing remaining in place for up to 4 days on a moderately exuding wound (1.1 g/cm²/24hrs). These calculations give an indication of the amount of fluid that is expected to be managed: actual values in use may be altered by environmental conditions, such as occlusion of the dressing or incorrect positioning of port.

References: 1. Smith+Nephew 2025. Calculation of the Theoretical Fluid Handling of PICO Dressings. Internal Report. RD/25/078.

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Advanced Wound Management
Smith & Nephew Medical Ltd
Croxley Park Building 5
Lakeside, Hatters Lane, Watford
Hertfordshire WD18 8YE, UK
T +44 (0)1482 225181
F +44 (0)1482 328326

www.smith-nephew.com

◊Trademark of Smith+Nephew
All Trademarks acknowledged
©October 2026 Smith+Nephew

AWM-AWD-49345 | GMC0952